



COMMODES

IMPERIAL COLLECTION
3-IN-1 DROP ARM COMMODE

Features

- Arms drop independently below seat level for safe lateral transfers
- Easy-to-release mechanism may be operated by patients lacking manual dexterity
- Arms automatically lock in raised position
- Extra wide platform seat is comfortable for long periods of sitting
- Welded steel construction and extra wide base provide added stability
- Cushioned armrests enhance user comfort
- Height readily adjusts from 19 ½" - 23 ½", in 1 inch increments
- Removable backbar facilitates use as a bedside commode, raised toilet seat or toilet safety frame
- Comes complete with commode pail, cover and splash shield
- Designed to meet the requirements of HCPCS code: E0168
- Maximum Weight Capacity 600 lb - EVENLY DISTRIBUTED
- Three (3) year limited warranty



Replacement Parts

		UOM
9290	Replacement Tips	4/bag

	overall Width		Seat Width		Seat Depth		minimum Seat Height		maximum Seat Height		Width between arms		Weight Capacity EVENLY DISTRIBUTED		UOM
	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	lb	Kg	
6438A	31	78.74	23	58.42	18	45.72	19.5	49.53	23.5	59.69	25	63.50	600	272.16	2/cs



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2008, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.