



COMMODOES

EVERYDAY COMMODE PAIL

Features

- Fits most Lumex[®] and other leading manufacturers' commodes
- 7 quart capacity
- Non-autoclavable
- Designed to meet the requirements of HCPCS code: E0167
- Made in U.S.A.



UOM

6691	7 qt commode pail, non-autoclavable	6/cs
PP600009	without cover	1 ea



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2022, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.