

COMMODOES

VERSAMODE™ DROP ARM

**Features**

- Provides transport, transfer, and commode functions all in one unit
- Arms drop independently below seat level for lateral transfers
- Cushioned armrests enhance user comfort
- Swing-away footrests can be extended or detached
- 5" rear swivel, locking casters facilitate maneuverability in tight spaces
- Padded seat cover
- Comes unassembled
- Designed to meet the requirements of HCPCS code: E0165
- Maximum Weight Capacity 300 lb - EVENLY DISTRIBUTED
- Three (3) year limited warranty



		overall Width		padded Seat Width		padded Seat Depth		commode Seat Width		commode Seat Depth		commode Seat Height		Width between arms		Weight Capacity EVENLY DISTRIBUTED		UOM
		inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	lb	Kg	
6810A	2 Rear Locking Casters	22	55.9	16	40.64	16	40.64	17.75	45.09	17.5	44.45	19.5	49.53	18.25	46.36	300	136.08	1 ea
6810B	4 Locking Casters	22	55.9	16	40.64	16	40.64	17.75	45.09	17.5	44.45	19.5	49.53	18.25	46.36	300	136.08	1 ea



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2012, GF Health Products, Inc. All Rights Reserved.  
Lumex is a registered trademark of GF Health Products, Inc.