



## COMMODOES

# PLATINUM COLLECTION 3-IN-1 STEEL COMMODOE

### Features

- Versatile 3-in-1 design can be used as a bedside commode, raised toilet seat or toilet safety frame
- Welded steel construction provides strength and stability
- Increased seat depth provides added stability and comfort
- Durable, white-powder coated steel frame is attractive and easy to clean
- Contoured waterfall armrests are designed to enhance user comfort and provide a good gripping surface
- Sturdy plastic seat snaps-on easily
- Seat height readily adjusts in 1 inch increments
- Comes complete with commode pail, cover and splash shield
- Designed to meet the requirements of HCPCS code: E0163
- Maximum Weight Capacity 300 lb - EVENLY DISTRIBUTED
- Three (3) year limited warranty



### Replacement Parts

		UOM
7103L	Replacement Legs	1 pr
7103S	Replacement Seat and Lid	1 ea
7103SG	Replacement Splash Guard	1 ea
RP20790-6	Replacement Pail and Lid	1 pr

7103	overall Width		overall Depth		Seat Width		Seat Depth		minimum Seat Height		maximum Seat Height		Width between arms		Weight Capacity EVENLY DISTRIBUTED		UOM
	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	lb	Kg	
	24.5	62.23	20	50.80	14	35.56	16.5	41.91	17	43.18	27	68.58	18.5	46.99	300	136.08	4/cs



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2009, GF Health Products, Inc. All Rights Reserved.  
Lumex is a registered trademark of GF Health Products, Inc.

graham  
field

GF0900159RevD22