



## COMMODES

# SILVER COLLECTION 3-IN-1 STEEL COMMODE

### Features

- Versatile 3-in-1 design can be used as a bedside commode, raised toilet seat or toilet safety frame
- Welded steel construction provides strength and stability
- Durable, grey-coated steel frame is attractive and easy to clean
- Contoured armrests enhance user comfort
- Backbar detaches to accommodate use as toilet safety frame
- Sturdy plastic seat snaps-on easily
- Seat height readily adjusts from 15 ½" to 23 ½", in 1 inch increments
- Comes complete with commode pail, cover and splash shield
- Designed to meet the requirements of HCPCS code: E0163
- Maximum Weight Capacity 300 lb - EVENLY DISTRIBUTED
- Five (5) year limited warranty



	overall Width		Seat Width		Seat Depth		minimum Seat Height		maximum Seat Height		Width between arms		Weight Capacity EVENLY DISTRIBUTED		UOM
	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	lb	Kg	
7103A-4	24	60.96	14	35.56	15	38.10	15.5	39.37	23.5	59.69	17.5	44.45	300	136.08	4/cs



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2012, GF Health Products, Inc. All Rights Reserved.  
Lumex is a registered trademark of GF Health Products, Inc.

graham  
field

GF-1200-026Rev1622