




COMMODES

3-IN-1 STEEL FOLDING COMMODE

Features

- Versatile 3-in-1 design can be used as a bedside commode, raised toilet seat or toilet safety frame
- Welded steel construction provides strength and stability
- Folds flat to less than 4" for convenient storage and travel
- Tool-free assembly
- Sturdy plastic seat snaps-on easily
- Comes complete with commode pail, cover and splash shield
- Maximum Weight Capacity 300 lb - EVENLY DISTRIBUTED
- Designed to meet the requirements of HCPCS code: E0163
- Five (5) year limited warranty



		overall Width		Seat Width		Seat Depth		minimum Seat Height		maximum Seat Height		Width between arms		Weight Capacity		UOM
		inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	lb	Kg	
7108A		24	60.96	14	35.56	16.5	41.91	17	43.18	23	58.42	20	50.80	300	136.08	4/cs
7108R-1	Retail Package	22	55.88	14	35.56	16.5	41.91	17	43.18	23	58.42	18	45.72	300	136.08	1 ea



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2007, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.

graham
field

GF0700098Rev1E2z