



COMMODES

BARIATRIC
STEEL FOLDING COMMODE

Features

- Versatile 3-in-1 design can be used as a bedside commode, raised toilet seat or toilet safety frame
- Welded steel construction provides strength and stability
- Folds flat for convenient storage and travel (to less than 7")
- Tool-free assembly
- Sturdy plastic seat snaps-on easily
- Comes complete with commode pail, cover and splash shield
- Designed to meet the requirements of HCPCS code: E0168
- Maximum Weight Capacity 650 lb - EVENLY DISTRIBUTED
- Three (3) year limited warranty



	overall Width		Seat Width		Seat Depth		minimum Seat Height		maximum Seat Height		Width between arms		Weight Capacity EVENLY DISTRIBUTED		UOM
	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	lb	Kg	
7109A-2	30	76.20	14	35.56	16.5	41.91	17.5	44.45	21.5	54.61	24.5	62.23	650	294.84	2/cs



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2010, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.