



COMMODES

IMPERIAL COLLECTION  
3-IN-1 STEEL COMMUNE

**Features**

- Extra wide platform seat is comfortable for long periods of sitting
- Commode opening is contoured and indented in the front to accommodate the needs of larger users
- Welded steel construction and extra wide base provide strength and stability
- Cushioned back bar enhances user comfort
- Seat height readily adjusts from 16 ½" to 23", in 1 inch increments
- Removable backbar facilitates use as a bedside commode, raised toilet seat or toilet safety frame
- Comes complete with commode pail, cover and splash shield
- Designed to meet the requirements of HCPCS code: E0168
- Maximum Weight Capacity 400 lb - EVENLY DISTRIBUTED
- Three (3) year limited warranty



**Replacement Parts**

		UOM
9140	Replacement Tips	4/bag

	overall Width		Seat Width		Seat Depth		minimum Seat Height		maximum Seat Height		Width between arms		Weight Capacity EVENLY DISTRIBUTED		UOM
	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	lb	Kg	
7446A-2	26	66.04	23	58.42	18	45.72	16.5	41.91	23.5	59.69	23	58.42	400	181.44	2/cs



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2012, GF Health Products, Inc. All Rights Reserved.  
Lumex is a registered trademark of GF Health Products, Inc.