



RAISED TOILET SEATS & SAFETY RAILS

ADJUSTABLE TOILET SAFETY RAIL

Features

- Designed to add support and enhance safety while toileting
- Contoured waterfall armrests are designed to add comfort and provide an excellent gripping surface
- Anodized aluminum frame is sturdy and lightweight
- Quickly secures to most toilets
- Arm width adjusts from 17" - 19"
- Easy tool-free assembly
- Designed to meet the requirements of HCPCS code: E0243
- Maximum Weight Capacity 250 lb - EVENLY DISTRIBUTED
- Six (6) month limited warranty



	overall Width		Width between arms		Weight Capacity EVENLY DISTRIBUTED		UOM
	inch	cm	inch	cm	lb	Kg	
6450R-2	22.5	57.15	17 - 19	43.18 - 48.26	250	113.40	2/cs



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2011, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.