



RAISED TOILET SEATS & SAFETY RAILS

ADJUSTABLE HEIGHT VERSAFRAME

**Features**

- Versaframe is designed to add support and enhance safety while toileting
- Contoured waterfall armrests are designed to add comfort and provide an excellent gripping surface
- Anodized aluminum frame is sturdy and lightweight
- Quickly secures to most toilets
- Arm height adjusts from 28" – 32"
- Arm width adjusts from 18" – 21"
- Easy tool-free assembly
- Designed to meet the requirements of HCPCS code: E0243
- Maximum Weight Capacity 250 lb - EVENLY DISTRIBUTED
- Six (6) month limited warranty



**Replacement Parts**

		UOM
LX9141	Replacement Tips	4/bag

		overall Width		minimum Arm Height		maximum Arm Height		minimum Arm Width		maximum Arm Width		Weight Capacity EVENLY DISTRIBUTED		UOM
		inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	lb	Kg	
6460A		24	60.96	28	71.12	32	81.28	18	45.72	21	53.34	250	113.40	4/cs
6460R	unassembled, retail package	24	60.96	28	71.12	32	81.28	18	45.72	21	53.34	250	113.40	3/cs
6465A-1	unassembled, brown box	24	60.96	28	71.12	32	81.28	18	45.72	21	53.34	250	113.40	1 ea



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2009, GF Health Products, Inc. All Rights Reserved.  
Lumex is a registered trademark of GF Health Products, Inc.