



RAISED TOILET SEATS & SAFETY RAILS

LOCKING RAISED TOILET SEAT

Features

- Designed for those who have difficulty sitting down or standing up from the toilet
- Adds 4 ½" to toilet seat height, reducing bending effort
- Front clamping mechanism locks securely to the toilet and prevents shifting of the seat
- Wide, contoured seating surface is comfortable
- White, heavy duty plastic construction is strong and durable
- Available in attractive retail package
- Designed to meet the requirements of HCPCS code: E0244
- Maximum Weight Capacity 300 lb - EVENLY DISTRIBUTED
- Six (6) month limited warranty



		overall Width		Seat Height		Seat Depth		Weight Capacity EVENLY DISTRIBUTED		UOM
		inch	cm	inch	cm	inch	cm	lb	Kg	
6486A		16.5	41.91	4.5	11.43	16.5	41.91	300	136.08	3/cs
6486R	Retail Packaging	16.5	41.91	4.5	11.43	16.5	41.91	300	136.08	3/cs



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2012, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.