



RAISED TOILET SEATS & SAFETY RAILS

LOCKING RAISED TOILET SEAT
WITH REMOVABLE ARMRESTS

Features

- Designed for those who have difficulty sitting down or standing up from the toilet
- Adds 4 1/2" to toilet seat height, reducing bending effort
- Front clamping mechanism locks securely to the toilet and prevents shifting of the seat
- Wide, contoured seating surface is comfortable
- White, heavy duty plastic construction is strong and durable
- Push buttons allow for easy removal of armrests
- Designed to meet the requirements of HCPCS code: E0244
- Maximum Weight Capacity 300 lb - EVENLY DISTRIBUTED
- Six (6) month limited warranty



	overall Width		Seat Height		Seat Depth		Width Between Arms		Weight Capacity EVENLY DISTRIBUTED		UOM
	inch	cm	inch	cm	inch	cm	inch	cm	lb	Kg	
6487RA	20.75	52.71	4.5	11.43	16.25	41.28	18.5	46.99	300	136.08	2/cs
6487RA-1	20.75	52.71	4.5	11.43	16.25	41.28	18.5	46.99	300	136.08	1 ea



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2007, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.