



RAISED TOILET SEATS & SAFETY RAILS

MULTI-POSITION  
OPEN PADDED RAISED TOILET

**Features**

- Unique open seat promotes patient independence and facilitates perineal cleaning and bowel training
- Front opening is 4 ¾" wide
- Seat can be rotated from front to either side or rear mounting for maximum access to perineal area
- Thick, heat-sealed padded seat is ideal for patients with circulatory problems or patients who need additional comfort
- Angle and height are adjustable to maximize comfort and reduce bending and patient effort required
- Comes complete with four heavy duty VersaGuard coated locking brackets
- Fits most toilets
- Optional extra-wide bracket (6499A) available for extra-wide toilet bowls
- Designed to meet the requirements of HCPCS code: E0244
- Maximum Weight Capacity 300 lb - EVENLY DISTRIBUTED
- Six (6) month limited warranty



	Seat Width		Seat Depth		minimum Seat Height		maximum Seat Height		Weight Capacity EVENLY DISTRIBUTED		UOM
	inch	cm	inch	cm	inch	cm	inch	cm	lb	Kg	
6497A	16	40.64	16	40.64	6	15.24	7.5	19.05	300	136.08	1 ea



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2012, GF Health Products, Inc. All Rights Reserved.  
Lumex is a registered trademark of GF Health Products, Inc.