



SHOWER TRANSPORT CHAIRS

SHOWER TRANSPORT CHAIR

Features

- Designed for safe transport to and from the shower
- Anodized, non-corrosive aluminum frame construction will not rust
- Comfortable nylon backrest removes for easy cleaning
- Cushioned armrests enhance user comfort
- 5" swivel casters with rear wheel safety locks facilitate maneuverability in tight spaces
- Rolls easily over most toilets
- Allows complete access to body for easy cleaning during bathing
- Maximum Weight Capacity
250 lb - EVENLY DISTRIBUTED
- Three (3) year limited warranty



	overall Height		overall Width		overall Depth		Seat Height		Seat Width		Seat Depth		Width between arms		Weight Capacity EVENLY DISTRIBUTED		UOM
	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	lb	Kg	
7910A-1	39	99	21.5	54.6	24	61	22	55.88	14	35.56	16.5	41.91	17.5	44.45	250	113.4	1 ea



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2012, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.