



SHOWER TRANSPORT CHAIRS

PVC KNOCKDOWN SHOWER CHAIR

Features

- Easy to assemble and breakdown for convenient transport
- Constructed of the highest quality health care grade PVC
- Deluxe elongated extra heavy-duty, contoured open front seat
- Dual reinforcement into the caster fittings adds additional support
- Quick drying, removable mesh back support is standard on all shower chairs
- Available in 3 anti-bacterial mesh or anti-bacterial solid vinyl colors
- Ergonomic curved arms provide support without sharp edges and include anti-slip handgrips for added safety
- Threaded stem casters for additional safety
- Clamp Joint Bracket Locking System is made of high grade impact plastic and stainless steel
- No sharp edges
- Will secure pipe and fittings 100%
- No tools required
- Extensively tested
- Patent pending on clamp joint bracket
- One (1) year limited warranty on PVC frame; six (6) month limited warranty on all other materials



Standard Mesh and Solid Vinyl Colors



Royal Blue



Forest Green



Maive

| | overall Height | | overall Width | | overall Depth | | Seat Height | | Seat Width | | Seat Depth | | Seat Clearance | | Width between arms | | Weight Capacity | | UOM |
|----------|----------------|--------|---------------|-------|---------------|-------|-------------|-------|------------|-------|------------|-------|----------------|------|--------------------|-------|-----------------|--------|------|
| | inch | cm | inch | cm | inch | cm | inch | cm | inch | cm | inch | cm | inch | cm | inch | cm | lb | Kg | |
| 89100-KD | 40.5 | 102.87 | 22.5 | 57.15 | 25.25 | 64.14 | 20 | 50.80 | 18 | 45.72 | 18 | 45.72 | 17.6 | 44.7 | 18 | 45.72 | 300 | 136.08 | 1 ea |
| 89200-KD | 40.5 | 102.87 | 26.5 | 67.31 | 27 | 68.58 | 20 | 50.80 | 22 | 55.88 | 18 | 45.72 | 17.6 | 44.7 | 22 | 55.88 | 375 | 170.10 | 1 ea |



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2008, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.

graham
field

GF0800210Rev/D22