



TRANSFER BENCHES

KNOCK DOWN TRANSFER BENCH

Features

- Provides assistance to those who have difficulty stepping over bathtub walls
- Anodized aluminum frame is lightweight, durable and rust-resistant
- Drainage holes in seat minimize water build-up
- Side support arm aids in lateral transfer
- Knock Down design is easy to assemble and saves space for transport and storage
- Seat height is adjustable in 1/2" increments
- Reversible back attachment to accommodate any bathroom
- Attractive platinum grey color
- Maximum Weight Capacity 350 lb - EVENLY DISTRIBUTED
- Five (5) year limited warranty



	overall Width		overall Depth		Seat Width		Seat Depth		minimum Seat Height		maximum Seat Height		minimum Clearance below frame		maximum Clearance below frame		Back Support height above seat		Weight Capacity EVENLY DISTRIBUTED		UOM
	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	lb	Kg	
7927KD-1	29.5	74.93	17	43.18	25.5	64.77	16	40.64	19.25	48.9	23.25	59.06	16.25	41.28	20.25	51.44	14.5	36.83	350	158.76	1 ea
7927KD-2	29.5	74.93	17	43.18	25.5	64.77	16	40.64	19.25	48.9	23.25	59.06	16.25	41.28	20.25	51.44	14.5	36.83	350	158.76	2/cs



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2015, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.

graham
field

GF-15-00059Rev1622