



## TRANSFER BENCHES

# IMPERIAL COLLECTION TRANSFER BENCH

### Features

- Provides assistance to those who have difficulty stepping over bathtub walls
- Seat has curtain cut-out and accessory tray
- Drainage holes in seat minimize water build-up
- Extra-deep seat provides added stability
- Anodized aluminum frame is lightweight, durable and rust-resistant
- Seat height is adjustable in ½" increments
- Reversible back attachment to accommodate any bathroom
- Attractive platinum grey color
- Comes assembled with tool-free back attachment
- Designed to meet the requirements of HCPCS code: E0248
- Maximum Weight Capacity 400 lb - EVENLY DISTRIBUTED
- Three (3) year limited warranty



	overall Width		overall Depth		Seat Width		Seat Depth		minimum Seat Height		maximum Seat Height		minimum Clearance below frame		maximum Clearance below frame		Back Support height above seat		Weight Capacity EVENLY DISTRIBUTED		UOM
	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	lb	Kg	
7929	32.5	82.55	23	58.42	28	71.12	16.75	42.55	17.5	44.45	21.5	54.61	15	38.1	19	48.26	13	33.02	400	181.44	2/cs
7929-1	32.5	82.55	23	58.42	28	71.12	16.75	42.55	17.5	44.45	21.5	54.61	15	38.1	19	48.26	13	33.02	400	181.44	1 ea



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2012, GF Health Products, Inc. All Rights Reserved.  
Lumex is a registered trademark of GF Health Products, Inc.

graham  
field

GF12000014RevB22