




TRANSFER BENCHES

PADDED KNOCK DOWN TRANSFER BENCH

Designed to provide assistance to those having difficulty stepping over bathtub walls.

Features

- Comfortable, cushioned seat and backrest are water-tight and easy to clean
- Anodized aluminum frame is lightweight, durable and rust-resistant
- Side support arm aids in lateral transfer
- Seat height is adjustable in 1" increments
- Reversible to accommodate any bathroom
- Knock-down design makes transporting easy
- Quick and easy tool-free assembly
- Maximum Weight Capacity 300 lb - EVENLY DISTRIBUTED
- Three (3) month limited warranty



		overall Width		overall Depth		Seat Width		Seat Depth		minimum Seat Height		maximum Seat Height		Commode diameter		Weight Capacity EVENLY DISTRIBUTED		UOM
		inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	lb	Kg	
7955KD-1	padded transfer bench	29	73.66	23	58.42	24	60.96	16	40.64	18	45.72	22	55.88	-	-	300	136.08	1 ea
7955KD-2	padded transfer bench	29	73.66	23	58.42	24	60.96	16	40.64	18	45.72	22	55.88	-	-	300	136.08	2/cs
7956KD-1	padded commode transfer bench	29	73.66	23	58.42	24	60.96	16	40.64	18	45.72	22	55.88	8	20.32	300	136.08	1 ea



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2018, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.

graham
field

GF-1800098RevC23