



HOMECARE BEDS

★★★ PATRIOT

Homecare Bed Series

SEMI-ELECTRIC HOMECARE BED



Slots in frame ensure proper rail placement



Locking pins secure frame to sleep deck



Designed and engineered by the same team who develop the GF Long Term Care beds.

Sleep surface and motors are **MADE IN THE USA!**

- The motor unit and bed design provide the ultimate in quiet, smooth operation
- Complete welded steel frame construction
- The head and foot sections can be raised or lowered simultaneously or independently with one hand control
- During a power failure the motor unit can lower both the head and foot section with just a single standard 9 volt battery
- Assembles easily with just five components
- Designed to meet the requirements of applicable HCPCS codes
- Total bed weight assembled: 166 lb (bed and motor only)
- Available with Fiberboard Decorative Walnut Bed-ends or High-impact Plastic Bed-ends
- Maximum patient weight 350 lbs
- 450 lb evenly distributed total weight capacity including patient, mattress, rails, and accessories
- Bed Height: Min. 15", Max. 24" (at deck level)
- Lift Range: Head - 60°, Knee - 30°
- Overall Bed Size: 87" L x 36" W
- No external pinch points!
- Slots in frame to ensure proper rail placement
- Locking pins secure bed ends to sleep deck Attached lanyards prevent loss while moving or storing
- Only homecare bed to offer a 4-inch extension kit that can extend either the head, foot, or both head and foot decking Bed length will be 84" with one kit, 88" with 2 kits
- The only homecare bed to offer optional corner mattress retainers

Grid Deck offers the following advantages:

- Much stronger frame
- Better mattress support
- No spring sag or bottoming out with heavier patients
- Easier to clean
- No connecting links or springs to lose



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com





HOMECARE BEDS



SEMI-ELECTRIC HOMECARE BED

| Fiberboard Decorative Walnut Bed-Ends | High Impact Plastic Bed-Ends | Description | UOM | Meets HCPCS |
|---------------------------------------|------------------------------|---|------|-------------|
| US0208 | US0208PL | Semi-Electric Bed Only | 1/ea | E0295 |
| US0208-RFPKG | US0208PL-RFPKG | Semi-Electric Bed with 1633-Rolled Foam Mattress and Full Chrome Rails | 1/ea | E0260 |
| US0208-RFPKGHR | US0208PL-RFPKGHR | Semi-Electric Bed with 1633-Rolled Foam Mattress and Half Chrome Rails | 1/ea | E0260 |
| US0208-RPKG | US0208PL-RPKG | Semi-Electric Bed with No Mattress and Full Chrome Rails | 1/ea | E0260 |
| US0208-RPKGHR | US0208PL-RPKGHR | Semi-Electric Bed with No Mattress and Half Chrome Rails | 1/ea | E0260 |
| US0218-RFPKG | US0218PL-RFPKG | Semi-Electric Bed with 1633- Rolled Foam Mattress and FDA Full Rails | 1/ea | E0260 |
| US0218-RFPKGHR | US0218PL-RFPKGHR | Semi-Electric Bed with 1633- Rolled Foam Mattress and FDA Half Rails | 1/ea | E0260 |
| US0218-RFPKGQR | US0218PL-RFPKGQR | Semi-Electric Bed with 1633- Rolled Foam Mattress and FDA Quarter Rails | 1/ea | E0260 |
| US0218-RPKG | US0218PL-RPKG | Semi-Electric Bed with No Mattress and FDA Full Rails | 1/ea | E0261 |
| US0218-RPKGHR | US0218PL-RPKGHR | Semi-Electric Bed with No Mattress and FDA Half Rails | 1/ea | E0261 |
| US0218-RPKGQR | US0218PL-RPKGQR | Semi-Electric Bed with No Mattress and FDA Quarter Rails | 1/ea | E0261 |

Warranty:

- Two year warranty for defects in workmanship and materials of mechanical components, frame, and electronics
- During the warranty period, defective items will be repaired or replaced at manufacturer's option. The warranty does not include any labor charges incurred in replacement part(s) installation or any associated freight or shipping charges to the manufacturer

Bed shipped in three cartons:

1. Foot section and motor.
2. Head section and Hi/Lo rod.
3. Bed end set.



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com

