



SUPPORT SURFACES

ADJUSTABLE BACK CUSHION

The E&J Adjustable Back Cushions allows for an individual fit for comfort and support.

Cushion Features

- 2" medium density foam
- 5 straps for maximum positioning options
- Wheelchair can still be folded by removing the back cushion

Cover Features

- Fluid-resistant stretch top
- Antimicrobial
- Flame retardant
- Easy to clean



Specifications

Item	UOM	Meets HCPCS	Overall Width		Overall Length		Overall Thickness		Made in USA	Limited Warranty in months
			inch	cm	inch	cm	inch	cm		
2055LS18	1/ea	E2611	18	45.7	17	43.2	2	5.1	■	12



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2013, GF Health Products, Inc. All Rights Reserved.
Everest & Jennings is a registered trademark of GF Health Products, Inc.

graham
field

GF1300007RevB23