



SUPPORT SURFACES

COMFORT CUSHION

The E&J Comfort Cushions provides pressure distribution and weight equalization.

Cushion Features

- 3" cushion
- Dual density foam for exceptional comfort and support base

Cover Features

- Vapour permeable (0% vapor transmission)
- Fluid resistant
- Infection Control is aided by antimicrobial properties; bacteria resistant, anti-fungal, and non-allergenic
- Anti-static and flame retardant
- Easy to remove cover with zipper handle for easy positioning
- Easy to clean and maintain



Specifications

Item	UOM	Meets HCPCS	Overall Width		Overall Length		Overall Thickness		Weight Capacity EVENLY DISTRIBUTED		Made in USA	Warranty in months
			inch	cm	inch	cm	inch	cm	lb	kg		
8100166	1/ea	E2611	16	40.6	16	40.6	3	7.6	300	136.1	■	12
8100168	1/ea	E2611	16	40.6	18	45.7	3	7.6	300	136.1	■	12
8100186	1/ea	E2611	18	45.7	16	40.6	3	7.6	300	136.1	■	12
8100188	1/ea	E2611	18	45.7	18	45.7	3	7.6	300	136.1	■	12
8100206	1/ea	E2611	20	50.8	16	40.6	3	7.6	300	136.1	■	12
8100208	1/ea	E2611	20	50.8	18	45.7	3	7.6	300	136.1	■	12



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2010, GF Health Products, Inc. All Rights Reserved.
Everest & Jennings is a registered trademark of GF Health Products, Inc.

graham
field

GF1000081Rev/D23