



WHEELCHAIR CUSHIONS



DURA-GEL®

SPP^{AT}

The Dura-Gel® center layer is surrounded by a comfortable, high-density foam base layer. The proprietary nonaqueous visco-elastic gel core will not leak and reduces tissue interface pressure at critical points.


Dura-Gel cushions are constructed with a fluid-resistant, removable and washable cover and are lightweight at under 5 pounds (2.27 Kg). The back securing strap, accommodates a wide range of body shapes, sizes, and weights. Packaged: 1 each



Features

- Antimicrobial
- Flame retardant
- Non-skid base
- Stretch fabric cover
- Visco-elastic foam topper
- Made in USA

★ See Chart for Item Numbers, Specifications and Warranty



SPP^{AT} Skin Protection & Positioning w/Anti-Thrust provides an extra level of positioning with an anti-thrust feature. A pre-ischial ridge design positions the user's back into the seat.

DURA-gel SPP^{AT} Hyperlink

Item	Meets HCPCS	Overall Width		Overall Length		Overall Thickness		Product ¹ Weight		Weight Capacity EVENLY DISTRIBUTED		Pressure Redistribution	Meets CTB117	Made in USA	Warranty in Months
		inch	cm	inch	cm	inch	cm	lb	kg	lb	kg				
8940186	E2607	18.00	45.72	16.00	40.64	3.00	7.62	4	1.82	300	136.08	■	■	■	24
8940208	E2607	20.00	50.80	20.00	50.80	3.00	7.62	4	1.82	300	136.08	■	■	■	24

■ Standard ■ Optional ■ Available — Not Applicable/Available Active Hyperlink
¹ approximate product weight — as size increases, so does product weight, however slightly



GF Health Products, Inc. is an Intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice. The most current and complete product information can be found on our website. www.grahamfield.com



© 2021, GF Health Products, Inc. All Rights Reserved. Everest & Jennings is a registered trademark of GF Health Products, Inc.