



WHEELCHAIR CUSHIONS




DURA-GEL®

SPIII

The Dura-Gel® center layer is surrounded by a comfortable, high-density foam base layer. The proprietary nonaqueous visco-elastic gel core will not leak and reduces tissue interface pressure at critical points.

Dura-Gel cushions are constructed with a fluid-resistant, removable and washable cover and are lightweight at under 5 pounds (2.27 Kg). The back securing strap, accommodates a wide range of body shapes, sizes, and weights.. Packaged: 1 each



Features

- Antimicrobial
- Flame retardant
- Non-skid base
- Stretch fabric cover
- Visco-elastic foam topper
- Made in USA

★ See Chart for Item Numbers, Specifications and Warranty

SPIII A great solution for pressure reduction/redistribution. A sloped, firm base foam has Dura-Gel units under the ischial tuberosity and upper thigh.

[DURA-gel SPIII Hyperlink](#)

Item	Meets HCPCS	Overall Width		Overall Length		Overall Thickness		Product ¹ Weight		Weight Capacity EVENLY DISTRIBUTED		Pressure Redistribution	Meets CTB117	Made in USA	Warranty in Months
		inch	cm	inch	cm	inch	cm	lb	kg	lb	kg				
8930186	E2603	18.00	45.72	16.00	40.64	3.00	7.62	5	2.27	350	158.76	■	■	■	18
8930208	E2603	20.00	50.80	20.00	50.80	3.00	7.62	5	2.27	350	158.76	■	■	■	18

■ Standard ■ Optional ■ Available — Not Applicable/Available [Active Hyperlink](#)
¹ approximate product weight — as size increases, so does product weight, however slightly



GF Health Products, Inc. is an Intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice. The most current and complete product information can be found on our website. www.grahamfield.com



© 2021, GF Health Products, Inc. All Rights Reserved. Everest & Jennings is a registered trademark of GF Health Products, Inc.