



WHEELCHAIR CUSHIONS




DURA-GEL®

BASE 3G

The Dura-Gel® center layer is surrounded by a comfortable, high-density foam base layer. The proprietary nonaqueous visco-elastic gel core will not leak and reduces tissue interface pressure at critical points.

Dura-Gel cushions are constructed with a fluid-resistant, removable and washable cover and are lightweight at under 5 pounds (2.27 Kg). The back securing strap, accommodates a wide range of body shapes, sizes, and weights.. Packaged: 1 each



Features

- Antimicrobial
 - Flame retardant
 - Non-skid base
 - Stretch fabric cover
 - Made in USA
- ★ See Chart for Item Numbers, Specifications and Warranty

BASE 3G The best general use cushion around and extremely lightweight. Supportive firm base foam and contouring top foam provide exceptional performance.

Item	Meets HCPCS	Overall Width		Overall Length		Overall Thickness		Product ¹ Weight		Weight Capacity EVENLY DISTRIBUTED		Pressure Redistribution	Meets CTB17	Made in USA	Warranty in Months
		inch	cm	inch	cm	inch	cm	lb	kg	lb	kg				
8565166	E2601	16.00	40.64	16.00	40.64	3.00	7.62	3	1.36	300	136.08	■	■	■	12
8565186	E2601	18.00	45.72	16.00	40.64	3.00	7.62	3	1.36	300	136.08	■	■	■	12
8565188	E2601	18.00	45.72	18.00	45.72	3.00	7.62	3	1.36	300	136.08	■	■	■	12
8565206	E2601	20.00	50.80	16.00	40.64	3.00	7.62	3	1.36	300	136.08	■	■	■	12
8565208	E2601	20.00	50.80	20.00	50.80	3.00	7.62	3	1.36	300	136.08	■	■	■	12
8565228	E2602	22.00	55.88	18.00	45.72	3.00	7.62	3	1.36	300	136.08	■	■	■	12
8565248	E2602	24.00	60.96	18.00	45.72	3.00	7.62	3	1.36	300	136.08	■	■	■	12

■ Standard ■ Optional ■ Available — Not Applicable/Available Active Hyperlink
¹ approximate product weight — as size increases, so does product weight, however slightly

DURA-gel BASE3G Hyperlink



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
 The most current and complete product information can be found on our website.
www.grahamfield.com



© 2021, GF Health Products, Inc. All Rights Reserved.
 Everest & Jennings is a registered trademark of GF Health Products, Inc.

graham
field

GF2100262 Rev B23