



ROLLATORS

WALKABOUT JUNIOR
FOUR-WHEEL ROLLATOR

This light-weight Walkabout Junior rollator is ideal for the pediatric user, has the same great features as the Walkabout Lite, weighs only 11.5 lbs, and is available in one attractive color. It folds quickly and easily into a compact unit for storage or transport. Its features include a padded seat, a backrest for extra comfort, ergonomic hand grips, and easy-to-operate hand brakes.

Features

- Epoxy-coated aluminum frame
- Locking loop brakes
- Removable wire storage basket
- Adjustable handle height
- 6" wheels for indoor or outdoor use
- Includes laminated color operating instructions
- Designed to meet the requirements of HCPCS codes: E0143 + E0156
- Three (3) year limited warranty



		Width inside handgrips		overall Width		overall Depth		Seat Width		Seat Depth		Seat Height	
		inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm
RJ4301R	Burgundy	12	30.48	20	50.80	24	60.96	12	30.48	12	30.48	18.5	46.99

		Wheel diameter		minimum Handle Height		maximum Handle Height		Product Weight		Weight Capacity EVENLY DISTRIBUTED		UOM
		inch	cm	inch	cm	inch	cm	lb	Kg	lb	Kg	
RJ4301R	Burgundy	6	15.24	27.5	69.85	32	81.28	11.5	5.22	300	136.08	1 ea



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2013, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.

Rollators are not intended to be propelled while seated.
Do NOT use this rollator as a wheelchair or transport device.