



## WHEELCHAIR CUSHIONS





# DURA-GEL®

BASE 2G

The Dura-Gel® center layer is surrounded by a comfortable, high-density foam base layer. The proprietary nonaqueous visco-elastic gel core will not leak and reduces tissue interface pressure at critical points.

Dura-Gel cushions are constructed with a fluid-resistant, removable and washable cover and are lightweight at under 5 pounds (2.27 Kg). The back securing strap, accommodates a wide range of body shapes, sizes, and weights. Packaged: 1 each



### Features

- Antimicrobial
- Flame retardant
- Non-skid base
- Stretch fabric cover
- Made in USA

★ See Chart for Item Numbers, Specifications and Warranty

**BASE 2G** Provides a durable yet comfortable solution for everyday seating needs. A layer of Dura-Gel and two densities of foam relieve pressure. Durable, 600 denier cover provides extra protection.

DURA-gel BASE2G Hyperlink

Item	Meets HCPCS	Overall Width		Overall Length		Overall Thickness		Product <sup>1</sup> Weight		Weight Capacity EVENLY DISTRIBUTED		Pressure Redistribution	Meets CTB117	Made in USA	Warranty in Months
		inch	cm	inch	cm	inch	cm	lb	kg	lb	kg				
8265166	E2601	16.00	40.64	16.00	40.64	2.00	5.08	3	1.36	250	113.34	■	■	■	12
8265186	E2601	18.00	45.72	16.00	40.64	2.00	5.08	3	1.36	250	113.34	■	■	■	12
8265188	E2601	18.00	45.72	18.00	45.72	2.00	5.08	3	1.36	250	113.34	■	■	■	12

■ Standard ■ Optional ■ Available — Not Applicable/Available Active Hyperlink  
<sup>1</sup> approximate product weight — as size increases, so does product weight, however slightly



GF Health Products, Inc. is an Intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice. The most current and complete product information can be found on our website. [www.grahamfield.com](http://www.grahamfield.com)



© 2021, GF Health Products, Inc. All Rights Reserved. Everest & Jennings is a registered trademark of GF Health Products, Inc.