



SUPPORT SURFACES

SIDE RAIL PADS

The Lumex Side Rail Pads are designed to help protect a patient from hard side rails by providing pressure relieving foam.

Features

- Designed to reduce gaps between hard bed rails and mattress
- Antimicrobial, flame and stain resistant vinyl cover
- Easy to clean
- Hook and loop fasteners attach easily to top rail
- For use with standard or swing-down bed rails



Specifications

Item	UOM	Overall Width		Overall Length		Standard Rail	Swing-Down Rail	Full Rail	Half Rail	3/4 Rail	Made in USA	Warranty in months
		inch	cm	inch	cm							
6013363	1/ea	14	35.6	36	91.4	■	—	—	■	—	■	12
6013364	1/ea	14	35.6	58	147.3	■	■	—	—	■	■	12
6013365	1/ea	17	43.2	30	76.2	—	■	—	■	—	■	12
6013663	1/ea	14	35.6	72	182.9	■	■	■	—	—	■	12