



## SUPPORT SURFACES


# ALTERNATING PRESSURE OVERLAY AND PUMP SYSTEM

### Features

- 2.5" (6.4 cm) bubble style overlay
- Heavy gauge, medical grade vinyl is durable, flame retardant, anti-bacterial
- Easy to clean
- Pump provides an automatic 5-minute alternating cycle time when it inflates and deflates.
- Available with fixed pressure or deluxe variable pressure pump
- Pump has built-in mounting brackets to easily mount to bed ends
- Pump has illuminated on/off switch



AQ1000



AQ2000



### Pump Specifications

	AQ1000	AQ2000
Power	120 Volt/60Hz, 1 Amp	120 Volt/60Hz, 1 Amp
Pressure	1.9 psi (13.1 kPa)	1.0 – 2.0 psi (6.9 – 13.8 kPa)
Air Output	4.5 liters/min	4.5 liters/min
Cycle Time	5 min	5 min
Power Cord Length	10 ft (3 m)	10 ft (3 m)
Dimensions	8.5" x 4.25" x 3.625" (21.6 x 10.8 x 9.2 cm)	9.75" x 5.125" x 3.875" (24.8 x 13 x 9.8 cm)
Weight	2.9 lb (1.3 kg)	2.5 lb (1.1 kg)

### Specifications

Item	Description	UOM	Meets HCPCS	Overall Width		Overall Length		Overall Height		Weight Capacity EVENLY DISTRIBUTED		Overlay in days	Pump in months
				inch	cm	inch	cm	inch	cm	lb	kg		
AQ1000	Standard pump and bubble-style pad, (120V)	1/ea	E0180	34	86.4	78	198.1	2.5	6.4	275	124.7	90	12
AQ1000-2	Standard pump only	1/ea	E0182	4.25	10.8	8.5	21.6	3.625	9.2	275	124.7	90	12
AQ1000-220V	Standard pump and bubble-style pad, (220V)	1/ea	E0180	34	86.4	78	198.1	2.5	6.4	275	124.7	90	12
AQ2000	Deluxe variable pressure pump and bubble-style pad	1/ea	E0181	34	86.4	78	198.1	2.5	6.4	275	124.7	90	12
AQ2000-2	Deluxe variable pressure pump only	1/ea	E0182	5.125	13.	9.75	24.8	3.875	9.8	275	124.7	90	12
AQ2000-1	Replacement pad for AQ1000 and AQ2000	1/ea	—	34	86.4	78	198.1	2.5	6.4	275	124.7	90	12



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2009, GF Health Products, Inc. All Rights Reserved.  
Lumex is a registered trademark of GF Health Products, Inc.

graham  
field

GF0900067/Rev/C23