




SUPPORT SURFACES

FLOOR MATS

The Lumex Floor Mats are designed to reduce injury by providing padding in the event a patient falls out of bed.

Features

- Designed for use with low beds
- Fold for convenient storage when not in use
- Available with beveled edge or Alzheimer bed style



BFP7224

Specifications

Item	UOM	Overall Width		Overall Length		Overall Thickness		Made in USA	Warranty in months	
		inch	cm	inch	cm	inch	cm			
BFP7224	Beveled Edge Floor Mat	1 ea	24	61	72	182.9	1.5	3.8	■	12
ECALZPADS	Alzheimer Bed Floor Mats	1 pr	30	76.2	72	182.9	2	5.1	■	12



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2012, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.