

EXAM ROOM SUPPLIES

STAINLESS STEEL
KICK BUCKET/STAND SET

Features

- Bucket pours easily without removal from stand
- Vinyl covered ring on stand prevents marring of walls and floors
- Smooth rolling 2" ball bearing conductive casters
- Capacity: 13 quarts
- Bucket dimensions: 11.8" D x 10" H
- Stainless steel



| | | bucket Diameter | | bucket Height | | bucket Volume | | caster Diameter | | UOM |
|--------|---------------------------|-----------------|----|---------------|------|---------------|------|-----------------|-----|------|
| | | inch | cm | inch | cm | qt | l | inch | cm | |
| 3267 | Kick Bucket and Stand Set | 11.81 | 30 | 10.04 | 23.5 | 13 | 12.3 | 2 | 5.1 | 1/ea |
| 3267 P | Utility Pail Only | 11.81 | 30 | 10.04 | 23.5 | 13 | 12.3 | — | — | 1/ea |



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2022, GF Health Products, Inc. All Rights Reserved.
Grafco is a registered trademark of GF Health Products, Inc.