



SUPPORT SURFACES

GEL-FOAM CUSHION

The E&J Gel-Foam Cushion provides pressure distribution and weight equalization.

Cushion Features

- 3" cushion
- High resiliency foam for support, comfort and durability
- Foam frame side walls provide support
- Gel channel bladders prevent bottoming out for maximum pressure relief
- Cover is removable and washable
- Retail packaged

Cover Features

- Vapour permeable (0% vapor transmission)
- Fluid resistant
- Infection control is aided by antimicrobial properties; bacteria resistant, anti-fungal, and non-allergenic
- Anti-static and flame retardant
- Easy to remove cover with zipper handle for easy positioning
- Easy to clean and maintain



Specifications

Item	UOM	Meets HCPCS	Overall Width		Overall Length		Overall Thickness		Weight Capacity EVENLY DISTRIBUTED		Warranty in months
			inch	cm	inch	cm	inch	cm	lb	kg	
81166	4 ea/cs	E2601	16	40.6	16	40.6	3	7.6	300	136.1	18
81186	4 ea/cs	E2601	18	45.7	16	40.6	3	7.6	300	136.1	18
81188	4 ea/cs	E2601	18	45.7	18	45.7	3	7.6	300	136.1	18
81216	4 ea/cs	E2601	20	50.8	16	40.6	3	7.6	300	136.1	18
81218	4 ea/cs	E2601	20	50.8	18	45.7	3	7.6	300	136.1	18
81220	4 ea/cs	E2601	20	50.8	20	50.8	3	7.6	300	136.1	18
81228	4 ea/cs	E2602	20	50.8	18	45.7	3	7.6	300	136.1	18
81248	4 ea/cs	E2602	24	61	18	45.7	3	7.6	300	136.1	18
81242	4 ea/cs	E2602	24	61	20	50.8	3	7.6	300	136.1	18



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2017, GF Health Products, Inc. All Rights Reserved.
Everest & Jennings is a registered trademark of GF Health Products, Inc.

graham
field

CF170003Rev1B23