


RECLINERS

REPLACEMENT PARTS


FR577RG CLINICAL CARE



23



24



25 Heat and Massage Feature – Control Module and Transformer Unit Showing All Connections Correctly Made



26



27

Item #	UOM	Description
5792A001G	1 ea	Table, Side
5792A002	4 ea/set	Hardware, Clevis Pin w Washer Side Table
57706J(color code)FR	1 ea	Arm Support, Upholstered
58735J(color code)FR	1 ea	Armrest, Upholstered
57702J(color code)FR	1 ea	Back, Upholstered
57702JH(color code)FR	1 ea	Back, Upholstered, Heat & Massage
57702JHO(color code)FR	1 ea	Back, Upholstered, Heat Only
57703J(color code)FR	1 ea	Footrest, Upholstered
57704J(color code)FR	1 ea	Legrest, Upholstered
57701J(color code)FR	1 ea	Seat, Upholstered
57701JH(color code)FR	1 ea	Seat, Upholstered, Heat & Massage
57701JHO(color code)FR	1 ea	Seat, Upholstered, Heat Only
56605J(color code)FR	1 ea	Side Panels, Upholstered
AP020004-8	8 ea/set	Hardware, Push Pin and Spring Side Panels
CSTRTRGN01	1 ea	Caster, Directional Green Tab Tente

Item #	UOM	Description
CSTRTRBLK01	1 ea	Caster, Swivel Locking Black Tab Tente
5772A216G	1 ea	Push Handle with Grip
9205A	1 ea	Push Handle with Grip Only for 5772A216G
5772A001	1 ea	Hardware, Cable Assembly
2740A020	1 pr	Hardware, Back Hinge
5772A300G	1 ea	Back Frame
5730F200	1 pr	Mechanism, Recliner
5740A023	1 pr	Hardware, Torque Tube
5792A003	10 ea/set	Hardware, Clevis Pin w Cotter Clip
5770A073	1 ea	Hardware, Gas Spring Bracket
EI000001	1 ea	Hardware, Gas Spring
LRHM0114	1 ea	Pendant Holster for Heat & Massage
HMOOG3CORD	1 ea	Heat & Massage Powercord Supply Gen 3
ECMCGF2	1 ea	Heat & Massage Control Box Gen 3
HMOOG3PEND	1 ea	Pendant, Heat & Massage Gen 3
HMOOG3PENDHO	1 ea	Pendant, Heat Only Gen 3
LRMP0001	1 ea	Paint, Touch-Up, Gray (not shown)



GF Health Products, Inc. is an ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

1 770.368.4700

Information contained herein is subject to change without notice. The most current and complete product information can be found on our website.

www.grahamfield.com



© 2018 GF Health Products, Inc. All Rights Reserved. Lumex is a trademark of GF Health Products, Inc. GF Health Products, Inc. is an ISO 13485 Certified Company.

graham field