

REPLACEMENT PARTS
FR566G DELUXE CLINICAL CARE



18

19

20

Heat and Massage Feature – Control Module and Transformer Unit Showing All Connections Correctly Made

Item #	UOM	Description
5792A001G	1 ea	Table, Side
5792A002	4 ea /set	Hardware, Clevis Pin w Washer, Side Table
5660C001G	1 ea	Trendelenburg Release Bar
PV860122	1 ea	Rubber Grip, Trendelenburg Bar
58735J(color code)FR	1 ea	Armrest, Upholstered
56606J(color code)FR	1 ea	Footrest, Upholstered
56604J(color code)FR	1 ea	Legrest, Upholstered
56601J(color code)FR	1 ea	Seat, Upholstered
56601JH(color code)FR	1 ea	Seat, Upholstered, Heat & Massage
56601JHO(color code)FR	1 ea	Seat, Upholstered, Heat Only
56605J(color code)FR	1 ea	Side Panels, Upholstered
AP020004-8	1 pr	Hardware, Push Pin and Spring, Side Panels
CSTRTRGRN01	1 ea	Caster, Directional Green Tab Tente
CSTRTRBLK01	1 ea	Caster, Swivel, Locking Black Tab Tente

Item #	UOM	Description
56602J(color code)FR	1 ea	Back, Upholstered
56602JH(color code)FR	1 ea	Back, Upholstered, Heat & Massage
56602JHO(color code)FR	1 ea	Back, Upholstered, Heat Only
PP460007	1 ea	Hardware, Back Attachment Thumb Screws
5772A216G	1 ea	Push Handle w Grip
9205A	1 ea	Grip Only, Push Handle
5660A200	1 ea	Mechanism, Recliner
LRHM0114	1 ea	Pendant Holster, Heat & Massage
HM00G3CORD	1 ea	Powercord Supply, Heat & Massage Gen 3
ECMCGF2	1 ea	Control Box, Heat & Massage Gen 3
HM00G3PEND	1 ea	Pendant, Heat & Massage Gen 3
HM00G3PENDHO	1 ea	Pendant, Heat Only Gen 3
LRMP0001		Paint, Touch-Up, Gray (not shown)



GF Health Products, Inc. is an ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice. The most current and complete product information can be found on our website.

www.grahamfield.com



© 2010, GF Health Products, Inc. All Rights Reserved. Lumex is a trademark of GF Health Products, Inc. GF Health Products, Inc. is an ISO 13485 Certified Company.

