



## WHEELCHAIRS



# PURETILT

The PureTilt® patented Tilt-in-Space Wheelchair, designed by an Occupational Therapist to provide safer and independent transitions to facilitate optimal comfort for a patient.

A manual tilt wheelchair that has a user-adjusting tilt control so that the wheelchair user can readjust their own posture throughout the day. This is vital for comfort and pressure relief.

The PureTilt comes standard with an ergonomic, contoured back that has adjustable height. Other standard features include anti-tippers for safety and a non removable padded leg rest that is length adjustable.



*folds for storage and vehicle transport*



# PURE TILT

## Features:

- User / caregiver-operated seat-tilt adjustment controlled by a built in lever that uses a gas spring piston
- Padded, hinged, non-detachable arms swing up and back to facilitate transfer
- Height-adjustable, ergonomic, contoured back with adjustable contoured headrest
- Back folds forward to facilitate storage and transport
- 22 inch rear wheels with aluminum hand rims for easy self- propulsion
- Two front casters that swivel for easy maneuvering
- Rear anti-tippers are a standard safety feature
- Angle-adjustable legrest assembly with the following features:
  - User-operated mechanism to adjust legrest angle
  - Footplate folds flat and locks in place to facilitate storage
  - Telescoping adjustable length

★ See Chart for Item Numbers, Specifications, and Warranty

## Upright postures help promote optimal:

- Food intake/digestion
- User comfortability during wheelchair propulsion
- Body positioning for comfort for increased out of bed tolerance

## Tilt in Space

- Provides user with more control over the position of their body in space to provide any pressure relief they may need
- Caregiver can also adjust positioning to help prevent discomfort when seated for extended periods
- Self-tilt promotes independence to participate in more daily/ leisurely activities



above shown with optional side pad set and foot pad accessories - each sold separately - see chart



padded, hinged, non-detachable arms shown in the up and back position at right

## Specifications: PT3000-18

|   | inch <sup>1</sup> | cm <sup>1</sup> |
|---|-------------------|-----------------|
| Overall Width:                                      | 28.50             | 72.39           |
| Overall Length: upright                             | 39.80             | 101.09          |
| Overall Length: reclined                            | 50.00             | 127.0           |
| Storage Length: back folded                         | 39.50             | 100.33          |
| Storage Height: back folded with seat               | 36.00             | 91.44           |
| Seat Width:   | 19.80             | 50.29           |
| Seat Depth:   | 21.00             | 53.34           |
| Seat Height: front of seat to floor - bottom pin    | 19.75             | 50.17           |
| Seat Height: front of seat to floor - top pin       | 21.75             | 55.25           |
| Arm Height: from seat - front                       | 5.80              | 14.73           |
| Foot Support Length: distance from seat - retracted | 13.00             | 33.02           |
| Foot Support Length: distance from seat - extended  | 17.50             | 44.45           |
| Front Caster Diameter:                              | 6.00              | 15.2            |
| Rear Wheel Diameter: with aluminum handrim          | 22.00             | 55.9            |

|   | degree of movement <sup>1,2</sup> |
|---|-----------------------------------|
| Seat Tilt: adjustment                   | 5 to 40                           |
| Leg to Seat: adjustment from horizontal | -97 to 15                         |

|                                       | lb  | kg <sup>1</sup> |
|---------------------------------------|-----|-----------------|
| Weight Capacity: EVENLY DISTRIBUTED   | 300 | 136.08          |
| Product Weight: includes all cushions | 95  | 43.09           |
| Shipping Weight: <sup>3</sup>         | 100 | 45.36           |

| Accessories:  | sold separately |
|---|-----------------|
| right & left Side Pads:<br>Item No. PT3000-SP       | ■               |
| larger leg pad plus Foot Pad:<br>Item No. PT3000-FP | ■               |

| Warranty:   | in months |
|---|-----------|
| Main, Seat & Back Frame:                          | 60        |
| Wheels, Casters, Legrest & Footplate:             | 12        |
| All Other Durable Components:<br>not listed above | 6         |

<sup>1</sup> ± 0.5 inches (± 1.27 cm) <sup>2</sup> ± 2° <sup>3</sup> estimated  
 † standard rounding method used to show metric values as whole numbers  
 ■ Standard ■ Optional — Does Not Apply



Standard Wheel Locks



Pin Adjustable Seat Height



Standard Adjustable Headrest



Optional Side Pad set (l&r)  
Item No. PT3000-SP  
sold separately



Forward-Folding Adjustable Back collapses for storage

GF Health Products, Inc. is an ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

• 1 770.368.4700

Information contained herein is subject to change without notice.

The most current and complete product information can be found on our website.  
[www.gendroninc.com](http://www.gendroninc.com) | [www.grahamfield.com](http://www.grahamfield.com)



© 2022, GF Health Products, Inc. All Rights Reserved.  
 Everest & Jennings is a trademark of GF Health Products, Inc.

**graham field**

GF21000073 Rev D23