



SUPPORT SURFACES



PATIENT POSITIONERS

The Lumex® Patient Positioners can be used for pressure relief, elevation, and specific patient positioning during procedures or recovery.

Features

- Made of 1.5 density foam for ultimate performance and comfort
- Latex free
- Meets California TB117 standard
- Non-sterile



P202




P204



P205



P301 - 14 INCH WIDTH
P302 - 16 INCH WIDTH
P303 - 18 INCH WIDTH



5004

Specifications		Overall Width		Overall Length		Overall Height		Shipping or Product Weight		UOM	Made in USA	Warranty in months
Item	Description	inch	cm	inch	cm	inch	cm	lb	kg			
P202	Arm Cradle	12.5	31.75	16.5	41.91	9	22.86	1.20	0.544	1/ea	■	6
P204	Heel & Elbow Protector ² Standard	6	15.24	7	17.78	2	5.08	0.35	0.159	1/pr	■	6
P205	Heel & Elbow Protector ² Boot-Shape	8	20.32	8.5	21.59	2	5.08	0.30	0.136	1/pr	■	6
P301	Abduction Pillow ² Small	14	35.56	18	45.72	5	12.70	1.17	0.531	1/ea	■	6
P302	Abduction Pillow ² Medium	16	40.64	22	55.88	6	15.24	1.42	0.644	1/ea	■	6
P303	Abduction Pillow ² Large	18	45.72	24	60.96	6	15.24	1.67	0.757	1/ea	■	6
5004	Body Aligner	12	30.48	24	60.96	6	15.24	1.13	0.513	1/ea	■	6

■ Standard ■ Optional ■ Available ■ Not Available ■ Accessory — Does Not Apply 🔗 Active Hyperlink
¹ approximate product weight — as seat width increases, so does product weight, however slightly ² with adjustable straps

[Product Webpage](#)



GF Health Products, Inc. is an Intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice. The most current and complete product information can be found on our website.

www.grahamfield.com



© 2012, GF Health Products, Inc. All Rights Reserved. Lumex is a registered trademark of GF Health Products, Inc.