



ROLLATORS

WALKABOUT IMPERIAL FOUR-WHEEL ROLLATOR

The Lumex Walkabout Four-Wheel Imperial (Bariatric) is ideal for larger individuals. Similar to the Walkabout Lite, this rollator offers a comfortable, lightweight and affordable solution to mobility. Designed with an aluminum frame, this lightweight rollator weighs 20 lbs and supports a maximum weight capacity of 450 lb - EVENLY DISTRIBUTED. The Walkabout folds quickly and easily into a compact unit for storage and transport.

Features

- Ergonomic hand grips
- Easy to operate locking and "quick stop" loop brakes (pull up to stop; push down to lock)
- Padded 20" wide seat for increased seating surface and added comfort
- Rear reflectors to increase visibility at night
- Removable padded backbar
- Adjustable handle height to accommodate broad range of user heights
- Durable, heavy duty 8" wheel for indoor and outdoor use on even terrain
- Pouch included that can be placed underneath seat for carrying of personal items
- Includes easy to read, retail-friendly user manual
- Designed to meet the requirements of HCPCS codes: E0149 & E1056
- Three (3) year limited warranty



		Width inside handgrips		overall Width		overall Depth		Seat Width		Seat Depth		Seat Height	
		inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm
RJ4400K	Straight Backbar, Black	21.5	54.61	28.7	72.90	26	66.04	19.5	49.53	13.8	35.05	21.8	55.37
RJ4400R	Straight Backbar, Burgundy	21.5	54.61	28.7	72.90	26	66.04	19.5	49.53	13.8	35.05	21.8	55.37
RJ4405B	Contoured Backbar, Royal Blue	21.5	54.61	28.7	72.90	26	66.04	19.5	49.53	13.8	35.05	21.8	55.37
RJ4405K	Contoured Backbar, Black	21.5	54.61	28.7	72.90	26	66.04	19.5	49.53	13.8	35.05	21.8	55.37
RJ4405R	Contoured Backbar, Burgundy	21.5	54.61	28.7	72.90	26	66.04	19.5	49.53	13.8	35.05	21.8	55.37

		Wheel diameter		minimum Handle Height		maximum Handle Height		Product Weight		Weight Capacity		UOM
		inch	cm	inch	cm	inch	cm	lb	Kg	lb	Kg	
RJ4400K	Straight Backbar, Black	8	20.32	34	86.36	39	99.06	20	9.07	450	204.12	1 ea
RJ4400R	Straight Backbar, Burgundy	8	20.32	34	86.36	39	99.06	20	9.07	450	204.12	1 ea
RJ4405B	Contoured Backbar, Royal Blue	8	20.32	34	86.36	39	99.06	20	9.07	450	204.12	1 ea
RJ4405K	Contoured Backbar, Black	8	20.32	34	86.36	39	99.06	20	9.07	450	204.12	1 ea
RJ4405R	Contoured Backbar, Burgundy	8	20.32	34	86.36	39	99.06	20	9.07	450	204.12	1 ea



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2009, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.

Rollators are not intended to be propelled while seated.
Do NOT use this rollator as a wheelchair or transport device.

graham
field

GF090006RevF23