



ROLLATORS

WALKABOUT WIDE FOUR-WHEEL ROLLATOR

18.5"
WIDE

The Lumex Walkabout Wide Four-Wheel Rollator offers a lightweight and affordable solution with now MORE added comfort. Designed with a wider aluminum frame, this lightweight rollator now offers a width inside the handgrips of 18.5" providing more room than that of traditional rollators without having to go to a bariatric model.

Weighing only 16 lbs, the rollator also supports up to 350 lb weight capacity. The Walkabout folds quickly and easily into a compact unit for storage and transport.

Features

- Ergonomic hand grips
- Easy to operate locking and "quick stop" loop brakes (pull up to stop; push down to lock)
- Wide padded seat for added comfort
- Removable padded straight backbar with easy release pull knobs
- Adjustable handle height to accommodate broad range of user heights
- 6" wheel for indoor and outdoor use
- Rear reflectors to increase visibility at night
- Basket included that can be placed underneath seat or in front for carrying of personal items
- Includes easy to read, retail-friendly user manual
- Designed to meet the requirements of HCPCS codes: E0143 & E1056
- Three (3) year limited warranty



		Width inside handgrips		overall Width		overall Depth		Seat Width		Seat Depth		Seat Height	
		inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm
RJ4318AQ	Aqua	18.5	46.99	27.5	69.85	27.5	69.85	15.5	39.37	13	33.02	22.5	57.15
RJ4318B	Royal Blue	18.5	46.99	27.5	69.85	27.5	69.85	15.5	39.37	13	33.02	22.5	57.15
RJ4318K	Black	18.5	46.99	27.5	69.85	27.5	69.85	15.5	39.37	13	33.02	22.5	57.15
RJ4318R	Burgundy	18.5	46.99	27.5	69.85	27.5	69.85	15.5	39.37	13	33.02	22.5	57.15

		Wheel diameter		minimum Handle Height		maximum Handle Height		Product Weight		Weight Capacity		UOM
		inch	cm	inch	cm	inch	cm	lb	Kg	lb	Kg	
RJ4318AQ	Aqua	6	15.24	32	81.28	37	93.98	16	7.26	350	158.76	1 ea
RJ4318B	Royal Blue	6	15.24	32	81.28	37	93.98	16	7.26	350	158.76	1 ea
RJ4318K	Black	6	15.24	32	81.28	37	93.98	16	7.26	350	158.76	1 ea
RJ4318R	Burgundy	6	15.24	32	81.28	37	93.98	16	7.26	350	158.76	1 ea



Rollators are not intended to be propelled while seated. DO NOT use this rollator as a wheelchair or transport device.



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2011, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.

graham
field

GF.1100039RevE23