




LIFT SLINGS

FULL-BODY MESH COMMODE SLING

Features

- Quick draining/drying mesh fabric offers full head and neck support
- For use with 4-point spreader option
- Designed to meet the requirements of HCPCS code: E0621
- Made in U.S.A.
- Not made with natural rubber latex
- Six (6) month limited warranty



Sling Size	Accent Color	body Width		body Length		Commode Opening to top of back		Weight Capacity EVENLY DISTRIBUTED		UOM	
		inch	cm	inch	cm	lb	Kg	lb	Kg		
FMC114	Medium	Red	42.75	108.59	50.75	128.91	34.75	88.27	450	204.1	1 ea
FMC115	Large	Green	44.75	113.67	54.25	137.8	38	96.52	450	204.1	1 ea
FMC116	X-Large	Blue	44.75	113.67	60.25	153.04	43	109.22	450	204.1	1 ea
FMC141	XX-Large	Yellow	42	106.68	51	129.54	41.5	105.41	600	272.2	1 ea



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2013, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.