





## LIFT SLINGS

# FULL-BODY MESH SLING

### Features

- Quick draining/drying mesh fabric offers full head and neck support
- For use with 4-point spreader option
- Designed to meet the requirements of HCPCS code: E0621
- Made in U.S.A.
- Not made with natural rubber latex
- Six (6) month limited warranty



| Sling Size | Accent Color | body Width |       | body Length |       | Weight Capacity<br>EVENLY DISTRIBUTED |     | UOM   |      |
|------------|--------------|------------|-------|-------------|-------|---------------------------------------|-----|-------|------|
|            |              | inch       | cm    | inch        | cm    | lb                                    | Kg  |       |      |
| FM110      | Medium       | Red        | 42.75 | 108.59      | 50.75 | 128.91                                | 450 | 204.1 | 1 ea |
| FM111      | Large        | Green      | 44.75 | 113.67      | 54.25 | 137.8                                 | 450 | 204.1 | 1 ea |
| FM140      | XX-Large     | Yellow     | 42    | 106.68      | 51    | 129.54                                | 600 | 272.2 | 1 ea |



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2013, GF Health Products, Inc. All Rights Reserved.  
Lumex is a registered trademark of GF Health Products, Inc.