



BARIATRIC PATIENT MOBILITY



# REGENCY WHEELCHAIR

## 6700 & 6800 SERIES FIXED-BACK & RECLINING WHEELCHAIRS

Designed for optimal bariatric patient transport as well as hospital and clinical use, Regency 6700 and 6800 wheelchairs are built with a dual square crossbrace heavy-duty frame and reinforced at all critical load-bearing points to support a wide variety of patients. With its multiple seat widths ranging from 22"-34", most with 18", 20" or 22" seat depths, and seat-height of 18", this heavy-duty wheelchair enables institutions and caregivers.

The reclining wheelchair series offers the flexibility to remain fully upright or recline close to a lying down position. This allows the patient the option to vary their position for greater comfort.

### Features

- Robust, steel reinforced frame
- Comfortable, easy-to-clean black nylon upholstery
- Removable arm styles with padded armrests
- Adjustable standard footrests or elevating legrest with durable black footplates
- Rear 24" x 1" wheels with black handrims and non-marking polyurethane tires
- 6700 models: Weight Capacity 700 lb - EVENLY DISTRIBUTED
- 6800 models: Weight Capacity 850 lb - EVENLY DISTRIBUTED
- Designed to meet the requirements of Medicare code: K0007 for fixed back models

### Reclining Model Only

- Pneumatic-assist reclining back to facilitate full recline positioning
- Removable full support headrest standard
- Adjustable length elevating legrest
- Rear anti-tippers standard
- Designed to meet the requirements of Medicare code: K0007/E1226 (Recliner Back)
  - + E0966 (Head Rest Extension)
  - + E0971 (Anti-tipping device)
  - + K0195 (Elevating Legrests)

### Warranty

- Five (5) year limited warranty



*top right: 6800 fixed-back model with full arms (shown) and elevating legrests (shown) — removable desk arm configuration available*

*right: 6800 fixed-back model with removable desk arms (shown) and adjustable swingaway footrest (shown) — other configurations available*







## Regency XL2000 Bariatric Reclining Wheelchair 6700 Series

Item	Meets HCPCS	Overall Chair Width		Overall Chair Depth		Overall Chair Height		Folded Chair Width		Seat Width		Seat Depth		Seat Height		Back Height		Wheel Diameter		Caster Diameter		Product <sup>1</sup> Weight w/o front rigging		Weight Capacity EVENLY DISTRIBUTED		Detach. Full Detach. Desk Adjust. Height Full <sup>3</sup>	Adjust. Swingaway Foot	Elevating Leg	Removable Head	Rear Anti-Tippers	Push Bar		
		inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	lb	kg	lb	kg	arms	rests								
67221830R	multi <sup>2</sup>	31	78.74	53	134.62	53	134.62	22	55.88	22	55.88	18	45.72	18	45.72	36	91.44	24	60.96	8	20.32	101	45.81	700	317.51	—	■	■	—	■	■	■	■
67221850R	multi <sup>2</sup>	31	78.74	53	134.62	53	134.62	22	55.88	22	55.88	18	45.72	18	45.72	36	91.44	24	60.96	8	20.32	101	45.81	700	317.51	■	—	■	—	■	■	■	■
67241830R	multi <sup>2</sup>	33	83.82	53	134.62	53	134.62	22	55.88	24	60.96	18	45.72	18	45.72	36	91.44	24	60.96	8	20.32	101	45.81	700	317.51	—	■	■	—	■	■	■	■
67241850R	multi <sup>2</sup>	33	83.82	53	134.62	53	134.62	22	55.88	24	60.96	18	45.72	18	45.72	36	91.44	24	60.96	8	20.32	101	45.81	700	317.51	■	—	■	—	■	■	■	■
67222030R	multi <sup>2</sup>	31	78.74	53	134.62	53	134.62	22	55.88	22	55.88	20	50.80	18	45.72	36	91.44	24	60.96	8	20.32	101	45.81	700	317.51	—	■	■	—	■	■	■	■
67222050R	multi <sup>2</sup>	31	78.74	53	134.62	53	134.62	22	55.88	22	55.88	20	50.80	18	45.72	36	91.44	24	60.96	8	20.32	101	45.81	700	317.51	■	—	■	—	■	■	■	■
67242030R	multi <sup>2</sup>	33	83.82	53	134.62	53	134.62	22	55.88	24	60.96	20	50.80	18	45.72	36	91.44	24	60.96	8	20.32	101	45.81	700	317.51	—	■	■	—	■	■	■	■
67242050R	multi <sup>2</sup>	33	83.82	53	134.62	53	134.62	22	55.88	24	60.96	20	50.80	18	45.72	36	91.44	24	60.96	8	20.32	101	45.81	700	317.51	■	—	■	—	■	■	■	■
67262030R	multi <sup>2</sup>	35	88.90	53	134.62	53	134.62	22	55.88	26	66.04	20	50.80	18	45.72	36	91.44	24	60.96	8	20.32	101	45.81	700	317.51	—	■	■	—	■	■	■	■
67262050R	multi <sup>2</sup>	35	88.90	53	134.62	53	134.62	22	55.88	26	66.04	20	50.80	18	45.72	36	91.44	24	60.96	8	20.32	101	45.81	700	317.51	■	—	■	—	■	■	■	■
67282030R	multi <sup>2</sup>	37	93.98	53	134.62	53	134.62	22	55.88	28	71.12	20	50.80	18	45.72	36	91.44	24	60.96	8	20.32	101	45.81	700	317.51	—	■	■	—	■	■	■	■
67282050R	multi <sup>2</sup>	37	93.98	53	134.62	53	134.62	22	55.88	28	71.12	20	50.80	18	45.72	36	91.44	24	60.96	8	20.32	101	45.81	700	317.51	■	—	■	—	■	■	■	■
67302030R	multi <sup>2</sup>	39	99.06	53	134.62	53	134.62	22	55.88	30	76.20	20	50.80	18	45.72	36	91.44	24	60.96	8	20.32	101	45.81	700	317.51	—	■	■	—	■	■	■	■
67302050R	multi <sup>2</sup>	39	99.06	53	134.62	53	134.62	22	55.88	30	76.20	20	50.80	18	45.72	36	91.44	24	60.96	8	20.32	101	45.81	700	317.51	■	—	■	—	■	■	■	■
67322030R	multi <sup>2</sup>	41	104.14	53	134.62	53	134.62	22	55.88	32	81.28	20	50.80	18	45.72	36	91.44	24	60.96	8	20.32	101	45.81	700	317.51	—	■	■	—	■	■	■	■
67322050R	multi <sup>2</sup>	41	104.14	53	134.62	53	134.62	22	55.88	32	81.28	20	50.80	18	45.72	36	91.44	24	60.96	8	20.32	101	45.81	700	317.51	■	—	■	—	■	■	■	■
67342030R	multi <sup>2</sup>	43	109.22	53	134.62	53	134.62	22	55.88	34	86.36	20	50.80	18	45.72	36	91.44	24	60.96	8	20.32	101	45.81	700	317.51	—	■	■	—	■	■	■	■
67342050R	multi <sup>2</sup>	43	109.22	53	134.62	53	134.62	22	55.88	34	86.36	20	50.80	18	45.72	36	91.44	24	60.96	8	20.32	101	45.81	700	317.51	■	—	■	—	■	■	■	■
67222230R	multi <sup>2</sup>	31	78.74	53	134.62	53	134.62	22	55.88	22	55.88	22	55.88	18	45.72	36	91.44	24	60.96	8	20.32	101	45.81	700	317.51	—	■	■	—	■	■	■	■
67222250R	multi <sup>2</sup>	31	78.74	53	134.62	53	134.62	22	55.88	22	55.88	22	55.88	18	45.72	36	91.44	24	60.96	8	20.32	101	45.81	700	317.51	■	—	■	—	■	■	■	■
67242230R	multi <sup>2</sup>	33	83.82	53	134.62	53	134.62	22	55.88	24	60.96	22	55.88	18	45.72	36	91.44	24	60.96	8	20.32	101	45.81	700	317.51	—	■	■	—	■	■	■	■
67242250R	multi <sup>2</sup>	33	83.82	53	134.62	53	134.62	22	55.88	24	60.96	22	55.88	18	45.72	36	91.44	24	60.96	8	20.32	101	45.81	700	317.51	■	—	■	—	■	■	■	■
67262230R	multi <sup>2</sup>	35	88.90	53	134.62	53	134.62	22	55.88	26	66.04	22	55.88	18	45.72	36	91.44	24	60.96	8	20.32	101	45.81	700	317.51	—	■	■	—	■	■	■	■
67262250R	multi <sup>2</sup>	35	88.90	53	134.62	53	134.62	22	55.88	26	66.04	22	55.88	18	45.72	36	91.44	24	60.96	8	20.32	101	45.81	700	317.51	■	—	■	—	■	■	■	■
67282230R	multi <sup>2</sup>	37	93.98	53	134.62	53	134.62	22	55.88	28	71.12	22	55.88	18	45.72	36	91.44	24	60.96	8	20.32	101	45.81	700	317.51	—	■	■	—	■	■	■	■
67282250R	multi <sup>2</sup>	37	93.98	53	134.62	53	134.62	22	55.88	28	71.12	22	55.88	18	45.72	36	91.44	24	60.96	8	20.32	101	45.81	700	317.51	■	—	■	—	■	■	■	■
67302230R	multi <sup>2</sup>	39	99.06	53	134.62	53	134.62	22	55.88	30	76.20	22	55.88	18	45.72	36	91.44	24	60.96	8	20.32	101	45.81	700	317.51	—	■	■	—	■	■	■	■
67302250R	multi <sup>2</sup>	39	99.06	53	134.62	53	134.62	22	55.88	30	76.20	22	55.88	18	45.72	36	91.44	24	60.96	8	20.32	101	45.81	700	317.51	■	—	■	—	■	■	■	■
67322230R	multi <sup>2</sup>	41	104.14	53	134.62	53	134.62	22	55.88	32	81.28	22	55.88	18	45.72	36	91.44	24	60.96	8	20.32	101	45.81	700	317.51	—	■	■	—	■	■	■	■
67322250R	multi <sup>2</sup>	41	104.14	53	134.62	53	134.62	22	55.88	32	81.28	22	55.88	18	45.72	36	91.44	24	60.96	8	20.32	101	45.81	700	317.51	■	—	■	—	■	■	■	■
67342230R	multi <sup>2</sup>	43	109.22	53	134.62	53	134.62	22	55.88	34	86.36	22	55.88	18	45.72	36	91.44	24	60.96	8	20.32	101	45.81	700	317.51	—	■	■	—	■	■	■	■
67342250R	multi <sup>2</sup>	43	109.22	53	134.62	53	134.62	22	55.88	34	86.36	22	55.88	18	45.72	36	91.44	24	60.96	8	20.32	101	45.81	700	317.51	■	—	■	—	■	■	■	■

■ Standard   ■ Optional   ■ Available   ■ Not Available   ■ Accessory   — Does Not Apply   🔗 Active Hyperlink

<sup>1</sup> approximate product weight — as seat width increases, so does product weight, however slightly  
<sup>2</sup> designed to meet the requirements of Medicare code: K0007/E1226 (Recliner Back) + E0966 (Head Rest Extension) + E0971 (Anti-tipping device) + K095 (Elevating Legrests)  
<sup>3</sup> adjustable height full arm accessory may be ordered separately at additional cost - item No. G10063

🔗 6700 Web Page



## Regency XL2000 Bariatric Fixed-Back Wheelchair

6800 Series

Item	Meets HCPCS	Overall Chair Width		Overall Chair Depth		Overall Chair Height		Folded Chair Width		Seat Width		Seat Depth		Seat Height		Back Height		Wheel Diameter		Caster Diameter		Product Weight w/o front rigging		Weight Capacity EVENLY DISTRIBUTED		Detach. Full arms	Detach. Desk	Adjust. Height Full <sup>3</sup>	Adjust. Swingaway Foot	Elevating Leg	Removable Head	Rear Anti-Tippers	Push Bar
		inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	lb	kg	lb	kg								
68262040	K0007	35	88.90	43	109.22	35	88.90	22	55.88	26	66.04	20	50.80	18	45.72	18	45.72	24	60.96	8	20.32	85	38.56	850	385.55	■	—	■	■	—	—	■	■
68282040	K0007	37	93.98	43	109.22	35	88.90	22	55.88	28	71.12	20	50.80	18	45.72	18	45.72	24	60.96	8	20.32	85	38.56	850	385.55	■	—	■	■	—	—	■	■
68302040	K0007	39	99.06	43	109.22	35	88.90	22	55.88	30	76.20	20	50.80	18	45.72	18	45.72	24	60.96	8	20.32	85	38.56	850	385.55	■	—	■	■	—	—	■	■
68322040	K0007	41	104.14	43	109.22	35	88.90	22	55.88	32	81.28	20	50.80	18	45.72	18	45.72	24	60.96	8	20.32	85	38.56	850	385.55	■	—	■	■	—	—	■	■
68342040	K0007	43	109.22	43	109.22	35	88.90	22	55.88	34	86.36	20	50.80	18	45.72	18	45.72	24	60.96	8	20.32	85	38.56	850	385.55	■	—	■	■	—	—	■	■
68262240	K0007	35	88.90	43	109.22	35	88.90	22	55.88	26	66.04	22	55.88	18	45.72	18	45.72	24	60.96	8	20.32	85	38.56	850	385.55	■	—	■	■	—	—	■	■
68282240	K0007	37	93.98	43	109.22	35	88.90	22	55.88	28	71.12	22	55.88	18	45.72	18	45.72	24	60.96	8	20.32	85	38.56	850	385.55	■	—	■	■	—	—	■	■
68302240	K0007	39	99.06	43	109.22	35	88.90	22	55.88	30	76.20	22	55.88	18	45.72	18	45.72	24	60.96	8	20.32	85	38.56	850	385.55	■	—	■	■	—	—	■	■
68322240	K0007	41	104.14	43	109.22	35	88.90	22	55.88	32	81.28	22	55.88	18	45.72	18	45.72	24	60.96	8	20.32	85	38.56	850	385.55	■	—	■	■	—	—	■	■
68342240	K0007	43	109.22	43	109.22	35	88.90	22	55.88	34	86.36	22	55.88	18	45.72	18	45.72	24	60.96	8	20.32	85	38.56	850	385.55	■	—	■	■	—	—	■	■
68262050	K0007	35	88.90	43	109.22	35	88.90	22	55.88	26	66.04	20	55.88	18	45.72	18	45.72	24	60.96	8	20.32	85	38.56	850	385.55	■	—	■	—	—	—	■	■
68282050	K0007	37	93.98	43	109.22	35	88.90	22	55.88	28	71.12	20	55.88	18	45.72	18	45.72	24	60.96	8	20.32	85	38.56	850	385.55	■	—	■	—	—	—	■	■
68302050	K0007	39	99.06	43	109.22	35	88.90	22	55.88	30	76.20	20	55.88	18	45.72	18	45.72	24	60.96	8	20.32	85	38.56	850	385.55	■	—	■	—	—	—	■	■
68322050	K0007	41	104.14	43	109.22	35	88.90	22	55.88	32	81.28	20	55.88	18	45.72	18	45.72	24	60.96	8	20.32	85	38.56	850	385.55	■	—	■	—	—	—	■	■
68342050	K0007	43	109.22	43	109.22	35	88.90	22	55.88	34	86.36	20	55.88	18	45.72	18	45.72	24	60.96	8	20.32	85	38.56	850	385.55	■	—	■	—	—	—	■	■
68262250	K0007	35	88.90	43	109.22	35	88.90	22	55.88	26	66.04	22	55.88	18	45.72	18	45.72	24	60.96	8	20.32	85	38.56	850	385.55	■	—	■	—	—	—	■	■
68282250	K0007	37	93.98	43	109.22	35	88.90	22	55.88	28	71.12	22	55.88	18	45.72	18	45.72	24	60.96	8	20.32	85	38.56	850	385.55	■	—	■	—	—	—	■	■
68302250	K0007	39	99.06	43	109.22	35	88.90	22	55.88	30	76.20	22	55.88	18	45.72	18	45.72	24	60.96	8	20.32	85	38.56	850	385.55	■	—	■	—	—	—	■	■
68322250	K0007	41	104.14	43	109.22	35	88.90	22	55.88	32	81.28	22	55.88	18	45.72	18	45.72	24	60.96	8	20.32	85	38.56	850	385.55	■	—	■	—	—	—	■	■
68342250	K0007	43	109.22	43	109.22	35	88.90	22	55.88	34	86.36	22	55.88	18	45.72	18	45.72	24	60.96	8	20.32	85	38.56	850	385.55	■	—	■	—	—	—	■	■

■ Standard 
 ■ Optional 
 ■ Available 
 ■ Not Available 
 ■ Accessory 
 — Does Not Apply 
 ⦿ Active Hyperlink

<sup>1</sup> approximate product weight — as seat width increases, so does product weight, however slightly  
<sup>3</sup> adjustable height full arm accessory may be ordered separately at additional cost - item No. G10063

[6800 Web Page](#)

## Regency XL2000 Bariatric Reclining Wheelchair 6800 Series

Item	Meets HCPCS	Overall Chair Width		Overall Chair Depth		Overall Chair Height		Folded Chair Width		Seat Width		Seat Depth		Seat Height		Back Height		Wheel Diameter		Caster Diameter		Product <sup>1</sup> Weight w/o front rigging		Weight Capacity EVENLY DISTRIBUTED		Detach. Full Arms	Detach. Desk	Adjust. Height Full <sup>3</sup>	Adjust. Swingaway Foot	Elevating Leg Rests	Removable Headrest	Rear Anti-Tippers	Push Bar
		inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	lb	kg	lb	kg								
68262050R	multi <sup>2</sup>	35	88.90	53	134.62	53	134.62	22	55.88	26	66.04	20	55.88	18	45.72	36	91.44	24	60.96	8	20.32	105	47.63	850	385.55	■	—	■	—	■	■	■	■
68282050R	multi <sup>2</sup>	37	93.98	53	134.62	53	134.62	22	55.88	28	71.12	20	55.88	18	45.72	36	91.44	24	60.96	8	20.32	105	47.63	850	385.55	■	—	■	—	■	■	■	■
68302050R	multi <sup>2</sup>	39	99.06	53	134.62	53	134.62	22	55.88	30	76.20	20	55.88	18	45.72	36	91.44	24	60.96	8	20.32	105	47.63	850	385.55	■	—	■	—	■	■	■	■
68322050R	multi <sup>2</sup>	41	104.14	53	134.62	53	134.62	22	55.88	32	81.28	20	55.88	18	45.72	36	91.44	24	60.96	8	20.32	105	47.63	850	385.55	■	—	■	—	■	■	■	■
68342050R	multi <sup>2</sup>	43	109.22	53	134.62	53	134.62	22	55.88	34	86.36	20	55.88	18	45.72	36	91.44	24	60.96	8	20.32	105	47.63	850	385.55	■	—	■	—	■	■	■	■
68262250R	multi <sup>2</sup>	35	88.90	53	134.62	53	134.62	22	55.88	26	66.04	22	55.88	18	45.72	36	91.44	24	60.96	8	20.32	105	47.63	850	385.55	■	—	■	—	■	■	■	■
68282250R	multi <sup>2</sup>	37	93.98	53	134.62	53	134.62	22	55.88	28	71.12	22	55.88	18	45.72	36	91.44	24	60.96	8	20.32	105	47.63	850	385.55	■	—	■	—	■	■	■	■
68302250R	multi <sup>2</sup>	39	99.06	53	134.62	53	134.62	22	55.88	30	76.20	22	55.88	18	45.72	36	91.44	24	60.96	8	20.32	105	47.63	850	385.55	■	—	■	—	■	■	■	■
68322250R	multi <sup>2</sup>	41	104.14	53	134.62	53	134.62	22	55.88	32	81.28	22	55.88	18	45.72	36	91.44	24	60.96	8	20.32	105	47.63	850	385.55	■	—	■	—	■	■	■	■
68342250R	multi <sup>2</sup>	43	109.22	53	134.62	53	134.62	22	55.88	34	86.36	22	55.88	18	45.72	36	91.44	24	60.96	8	20.32	105	47.63	850	385.55	■	—	■	—	■	■	■	■

■ Standard 
 ■ Optional 
 ■ Available 
 ■ Not Available 
 ■ Accessory 
 — Does Not Apply 
 🔗 Active Hyperlink

<sup>1</sup> approximate product weight — as seat width increases, so does product weight, however slightly  
<sup>2</sup> designed to meet the requirements of Medicare code: K0007/E1226 (Recliner Back) + E0966 (Head Rest Extension) + E0971 (Anti-tipping device) + K095 (Elevating Legrests)  
<sup>3</sup> adjustable height full arm accessory may be ordered separately at additional cost - item No. G10063

[6800 Web Page](#)



GF Health Products, Inc. is an ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

**+1 770.368.4700**

Information contained herein is subject to change without notice. The most current and complete product information can be found on our website.  
[www.gendroninc.com](http://www.gendroninc.com) | [www.grahamfield.com](http://www.grahamfield.com)



© 2021, GF Health Products, Inc. All Rights Reserved. Gendron is a registered trademark of GF Health Products, Inc.



GF2100003 Rev B23