



SUPPORT SURFACES

SELECT FOAM MATTRESSES

LS100 Series

- High Density Foam Mattress
- Convoluted Topper for Excellent Airflow
- Durable Nylon/Vinyl Cover
- Compress rolled for cost effective storage and freight
- ★ See Chart below for Links, Item Numbers, Specifications, Options and Warranty

LS150 Series

- High Density Foam Mattress
- Firm Side Rails and Soft Heel Section
- Long Term Care Replacement Mattress
- Stretch Urethane/Vinyl Cover
- Compress rolled for cost effective storage and freight
- Available with Raised Perimeter Bolster (LS150-35P and LS150-42P)
- ★ See Chart below for Links, Item Numbers, Specifications, Options and Warranty



Raised Perimeter Bolster
for patients likely to fall out of or exit the bed unattended

Firm Side Rails
aids in fall prevention during ingress / egress

Convoluted Topper
creates excellent airflow

Soft Heel Section
for patients with sensitive feet

Specifications

	Meets HCPCS	Width		Length		Height		Weight Capacity		Durable Nylon/Vinyl cover	Stretch Urethane/Vinyl	High Density Foam	Convoluted Topper	Firm Side Rails	Soft Heel Section	Raised Perimeter	Pressure Redistribution	Meets 16 CFR 1633	Mattress Warranty in Months	Cover Warranty in Months
		inch	cm	inch	cm	inch	cm	lb	kg											
LS100-35	E0184	35.00	88.90	80.00	203.20	6.00	15.24	350	158.76	■	—	■	■	—	—	—	■	■	60	12
LS100-42	E0184	42.00	106.68	80.00	203.20	6.00	15.24	350	158.76	■	—	■	■	—	—	—	■	■	60	12
LS150-35	E0184	35.00	88.90	80.00	203.20	6.00	15.24	400	181.44	—	■	■	—	■	■	—	■	■	60	12
LS150-42	E0184	42.00	106.68	80.00	203.20	6.00	15.24	400	181.44	—	■	■	—	■	■	—	■	■	60	12
LS150-35P	E0184	35.00	88.90	80.00	203.20	6.00	15.24	400	181.44	—	■	■	—	■	■	—	■	■	60	12
LS150-42P	E0184	42.00	106.68	80.00	203.20	6.00	15.24	400	181.44	—	■	■	—	■	■	—	■	■	60	12

■ Denotes Standard Feature

— Denotes Feature Not Available / Does Not Apply



Shipping Box Dimensions:
40 x 14 x 14 inches
(101.60 x 35.56 x 35.56 cm)



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2018, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.