



BARIATRIC PATIENT TRANSPORT



HEAVY DUTY TRANSPORT STRETCHER

The Gendron Heavy Duty Transport Stretcher is designed for treatment and transport of larger patients. Proudly made in the USA, this Heavy Duty Stretcher offers a weight capacity of 1,000 lbs evenly distributed. Gendron stretchers are designed with features that provide effortless and reliable patient handling.

**Features**

- Hydraulic Fowler backrest adjusts from 0-80°
- Steel powder-coat painted construction
- Two-position vertical side rails
- Heavy duty 8" casters with steer and brake functions
- 4" high-density pad with nylon cover
- Four IV sockets — one at each corner of stretcher



**Dimensions**

Stretcher Model	Patient Weight Capacity		overall Width		overall Length		Patient Surface (L x W)		Stretcher Top Height		Side Rail Length		Side Rail Height above Stretcher Top		Caster Size		Fowler Angle	Stretcher Shipping Weight	
	lb	kg	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	degrees	lb	kg
	EVENLY DISTRIBUTED																		
G1190	1,000	453.6	40	101.6	80	203.2	79 x 39	200.7 x 99.1	29.5	74.9	57.5	146.1	10.5	26.7	8	20.3	≤80	330	149.7



GF Health Products, Inc. is an Intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice. The most current and complete product information can be found on our website.

[www.gendroninc.com](http://www.gendroninc.com) | [www.grahamfield.com](http://www.grahamfield.com)



© 2020, GF Health Products, Inc. All Rights Reserved. Gendron is a registered trademark of GF Health Products, Inc.