

MAXI REST BARIATRIC BED

MAXI REST 4054 SERIES

Features

- Robust steel two-piece split-frame construction for ease of delivery and transport
- Bed deck width is expandable from 39" to 48" to 54"
- Bed deck length is adjustable from 80" to 88"
- 5" Deluxe Caster Set includes 2 Total-Lock and 2 Steer-Lock for ease of maneuverability
- Standard 3-Function Positioning hand control
- Optional 5-Function Positioning hand control (Trendelenburg / Reverse Trendelenburg and Chair positioning)
- Synchronized Hi/Lo Motors ensure bed remains level
- Battery backup standard for dependability
- Optional Intuitive side rail system for both head and foot sections
- Optional laminate wood grain or plastic molded head and foot board options
- Optional factory-installed weight scale system
- Safe working load: 1,000 lb (453.6 kg)*



Dimensions & Specifications

Bed Model	Patient Weight Capacity EVENLY DISTRIBUTED		Height Range		overall Length		Sleep Deck Width Range		Number of Functions	Rolling Position	Number of Motors	Built-in Scale	Emergency Battery Backup
	lb	kg	inch	cm	inch	cm	inch	cm	functions	height	motors	scale	battery
4054	1,000	453.6	13.5-27.5	34.3-69.9	80 or 88	203.2 or 223.5	39-48 or 54	99.1-121.9 or 137.2	3 (std) or 5 (optional)	any height	four	optional	yes

* Refer to user manual

MAXI REST 4054 SERIES




Maxi Rest Options and Accessories

available at additional cost

- GL928000 38 inch (96.52 cm) Laminate Wood Grain Head and Foot Board Set
- GL929000 39 inch (99.06 cm) HDPE Plastic Head and Foot Board Set
- GL930000 48 inch (121.92 cm) HDPE Plastic Head and Foot Board Set
- GL927100 Deluxe Plastic Rail Set (LH/RH) [4054 ONLY]
- GL927020 Standard Painted Foot End Rails w/ 39, 42 & 48 inch Pull-outs
- GL927000 Standard Painted Head End Rails
- GL93100 80 inch (203.20 cm) Frame Mount Full Length Trapeze [FITS 80 INCH ONLY]
- GL93180 88 inch (223.52) Frame Mount Full Length Trapeze [FITS 88 INCH ONLY]
- GL93200 5-Function Hand Control (Trend / Rev Trend / Chair)
- GL93300 Full Width Wall Bumper
- GL93600L Egress (Assist) Bar and Left Hand Mounting Pullout
- GL93600R Egress (Assist) Bar and Right Hand Mounting Pullout
- GL94000 8 inch (20.32 cm) Bolster Mattress and Extension Pan Kit
- GS4054DW Bed with Factory Installed Scale System


GL929000
39 INCH HDPE PLASTIC HEAD AND FOOT BOARD SET



GL927000
STANDARD PAINTED RAIL SET SHOWN WITH HAND CONTROL



GL927100
DELUXE PLASTIC RAIL SET SHOWN WITH HAND CONTROL



Warranty

- Main Frame and welds:..... Five years
- Control Box and Actuator Motors:..... Two years
- Hand Control Pendant, Staff Control, and Cabling:..... Two years
- Head Board and Foot Board:..... One year
- All other durable components not listed above:..... Two years

GL93200
5-FUNCTION TREND/REV
TREND HAND CONTROL



GL93100
FRAME MOUNT FULL LENGTH TRAPEZE
80 inch (203.20 cm) version shown



GS4054DW
BED WITH FACTORY INSTALLED SCALE SYSTEM



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.gendroninc.com | www.grahamfield.com



© 2020, GF Health Products, Inc. All Rights Reserved.
Gendron is a registered trademark of GF Health Products, Inc.

graham
field

GF200023Rev B23