



PAIN MANAGEMENT

TENS is a device that uses electrical current for stimulating nerves to help relieve pain. Usually used as an alternative to medications.

T.E.N.S. (Transcutaneous Electrical Nerve Stimulation) is indicated under a physician's prescription for the symptomatic relief and management of chronic (long-term) pain and the treatment of post-operative or post-traumatic pain



GF-3

Analog TENS

Prescription Required



- Feature:** Dual Channel
- Modes:** 3 - Standard, Burst, Modulation
- Includes:** TENS Unit, soft carry case, instructions, four electrodes and lead wires
- Uses:** 9 Volt Battery (NOT INCLUDED)

GF-3T

Analog TENS with Timer

Prescription Required



- Features:** Dual Channel, Timer
- Modes:** 3 - Standard, Burst, Modulation
- Includes:** TENS Unit, soft carry case, instructions, four electrodes and lead wires
- Uses:** 9 Volt Battery (NOT INCLUDED)

GF-TX5

Digital TENS

Prescription Required



- Features:** Dual Channel, Timer, LCD Screen
- Modes:** 3 - Standard, Burst, Modulation
- Includes:** TENS Unit, soft carry case, instructions, four electrodes and lead wires
- Uses:** 9 Volt Battery (NOT INCLUDED)



GF Health Products, Inc. is an ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice. The most current and complete product information can be found on our website.

www.grahamfield.com



© 2021, GF Health Products, Inc. All Rights Reserved.