

EXTENDED CARE BEDS



We proudly make all of our beds in the USA to the highest safety standards.<sup>6</sup>

- **FDA Recognized Standard:** ANSI/AAMI STD ES60601-1
- **Health Canada Recognized Standard:** CAN/CSA C22.2 No. 60601-1 (IEC 60601-1:2012-Edition 3.1)
- **Includes International Standards:** IEC 60601-1, IEC 60601-1-2, IEC 60601-2-52

## All Basic American Extended Care Beds Are

- Certified transport device - passes threshold clearance at ALL heights - ideal for evacuation
- Silver Solutions powder-coating provides antimicrobial surface protection and deters the growth of odor-causing bacteria
- No holes on top surface of tubes that could harbor fluids/bacteria
- Grounded bed with right-angle hospital grade plug for safety
- Adjustable Wallsaver protects the wall with smallest footprint
- Patented Hi-Lo mechanism adds life to the motors
- Sturdy, no-sag, grid sleep surface with cross braces
- Bed provides 360 mobility for cleaning, etc.
- Proudly Made in the USA
- Assist devices available
- Tool-less assembly



### Adjustable Wallsaver

- Maximum wall protection
- Adjustable clearance lengths
- Standard on all beds

Optional  
Zenith 7200  
Pendant Shown



### Mattress Retainers

- Adjustable for multiple mattress lengths
- Standard on all beds

## Zenith & Matrix Series Bed Standard Features

- Advanced Positioning System (APS) includes One-Button Comfort/Cardiac Chair Position (optional on Matrix 7200).
- Straight-Line Hi-Lo Travel with Patent-Awarded Design Technology
- 'Plug and Play' Digital Electronics enable future feature additions
- Smart Technology
- Industry-leading, energy-saving electronics consume up to 66% less power than transformer electronics when idle
- Sealed ball bearings joints provide continuous, maintenance-free, silent operation
- Intuitive two-pedal locking system - centrally located, easily visible and accessible at all heights and widths
- Tool-less assembly

## Zenith & Matrix Series Bed Configurations

- Embedded staff control with safety-inspired, two-stage lockouts enabling caregivers to lock-out features to the resident
- Standard or Advanced Positioning System (APS)
- APS Staff Controls are optional. (APS Staff Control includes One-Button Comfort/Cardiac Chair and Trendelenburg/Reverse Trendelenburg Positions)
- Hand Pendant to control Trendelenburg/Reverse Trendelenburg (*sold separately*)
- Accessories and some options are sold separately



### Intuitive Two-Pedal Locking System

- Centrally positioned and easily visible from a distance ensures the caregiver that bed is locked at a glance.
- Bed rolls at any height increasing caregiver safety
- Readily accessible at ALL bed widths and heights (even lowest)
- Select Option at time of purchase

# Zenith

## Zenith Series Standard Features

- With US Patent awarded time-saver Hi-Lo system
- Tall, Tubular Securely Mounted Head and Footboards create a built-in push bar and reduce breakage

## Zenith 9200 SERIES

9200 WITH APS SHOWN WITH  
Slide-W-I-D-E® deck option



## Zenith 7200 SERIES

7200 & 7200 WITH APS



## OPTIONAL WITH THE ZENITH SERIES *Slide-W-I-D-E*® EXPANDABLE/RETRACTABLE SLEEP DECK



Quick-release  
ball-detent pins  
attached to beds  
with lanyards



### Wider sleep surfaces reduce falls by residents accustomed to sleeping on common consumer products.<sup>2</sup>

- Patented design allows for quick adjustment, ideal for moving through doorways to assist in emergency evacuations
- Expandable/retractable deck quickly and easily adjusts to 35, 39, or 42 inches with two quick-release ball-detent pins on each side.
- Wide head and footboards available in 42 inch width
- All pins are attached to the bed with lanyards, eliminating lost parts
- Rails and assist devices are also quick-release without tools

<sup>2</sup> Fragala G, Perry B, Fragala M. Examining bed width as a contributor to risk of falls from bed in long term care. *Annals of Long-Term Care: Clinical Care and Aging*. 2012; 20(6):35-38



**Matrix 6200 APS** SHOWN  
with intuitive pedal-lock option  
Raised Head Position

# Matrix

## Matrix Series Standard Features

- Advanced Positioning System (APS) includes One-Button Comfort/Cardiac Chair Position (optional on Zenith 7200).
- Straight-Line Hi-Lo Travel with Patent-Awarded Design Technology

## Matrix 6200 SERIES 6200 WITH APS

## Advanced Positioning System (APS)

Advanced positioning allows the bed to go into numerous positions

- Positively influencing clinically driven outcomes through multiple positioning options
- Accommodating post-surgical rehab positioning to facilitate resident comfort
- Reduces caregiver strain and potential injuries versus manually repositioning residents
- Built-In Hi-Lo Lockout Feature (on APS Hand Controls)



**Zenith 9200 APS** SHOWN  
with Advanced Positioning System, Slide-W-I-D-E®  
deck option, and intuitive pedal-lock option  
Comfort/Cardiac Chair Position



**Liberty** SHOWN  
with intuitive pedal-lock option  
Supine Position



## Liberty Bed Standard Features


- Desired features at affordable price point
- Four Locking casters standard

## Liberty Bed Available with


- Maximize security with optional one-step red/green pedal lock (*as shown*)
- Optional 4-inch Total-lock casters increase under-bed clearance

## Overbed Table Options






## Counter-Rotating Assist Device




Provide a sturdy and secure hand hold to assist residents in and out of bed yet easily swivels out of the way when desired.

- Meet FDA Entrapment Guidelines
- Patented Counter-Rotating Assist Devices
- Removable without tools
- Move out of the way for lateral transfers
- Pivot away from egress area
- Provide toe clearance in lowest position to prevent staff injury
- Hand pendant mounts facing inward or outward and away from hand grab area



## Pivot Assist Bars



## Fixed Assist Devices

Provide a sturdy and secure hand hold to assist residents in and out of bed.

- Move out of the way for lateral transfers
- Tool-free assembly and removal



## Bed Cart

- Ideal for ease of storage.
- Bed cart allows beds to fit through narrow doorways.
- Designed for One-Person Mounting and maneuverability
- Works with all beds



## USB Power Supply

- Charges most devices
- Works with all beds



## Side Wall Bumper

- Maximum wall protection
- Works with all Zenith and Matrix Beds


## Advanced Positioning Pendant

- Upgrade from Standard to APS on Zenith 7200
- One-Button Comfort/Cardiac Chair Position
- Hi-Lo Lockout Feature



*Custom Finishes available upon request. Additional cost may incur.*





## Board Mounted IV Socket

- Mounts with 2 screws
- Easy to install
- IV pole sold separately
- Works with all beds




## Plastic Pan Deck Kit

- Fits 6200, 7200, 9200 Series



## Edema Foot Ratchet Kit

- Elevated foot rest.
- Standard on Zenith 9200
- Optional on Matrix 6200 and Zenith 7200



## Surface Length Extension Kit

- Extends standard length bed frame from 84" to 88". Complete with solid pan extension, mattress guard, and installation hardware.
- Works with all beds
- Easy to install




## UnderBed Light

- Reduces fall risks
- Eliminates bright light intrusion
- Multi-level light intensity
- Works with all beds



## Pin-on Wide Kit

- Fits 35 inch width 6200, 7200, 9200
- Available in 39" or 42" width kits
- Anytime add-on



## Embedded Staff Control

- Standard Functionality and Lockouts, or Advanced Positioning (APS).
- Angled-Up Display for Improved Visibility
- Two-Stage lockouts deter unwanted changes
- Mix or Match with Plug & Play design
- Upgrade anytime

## Trapeze and Adapter

- Pivots out of the way when not in use
- Works with all Zenith and Matrix Beds
- Rigid Design

Accessories

Kensington  
Maple

Solar Oak

Fonthill Pear

Wild Cherry

Windsor  
Mahogany

Shaker Cherry

Talas Cherry

Aria

Cocobala

## Bed Specifications:

	Zenith 9200	Zenith 7200	Matrix 6200	Liberty
Adjustable Patient Surface Length:	80-84 inch	80-84 inch	80-84 inch	76-80 inch
Welded Sleep Deck Width:	35 inch	35 inch	35 or 42 inch	35 inch
Height Travel Range:	7-30 inch	7.95-30 inch	8.95-27 inch	9.95-24.4 inch
Safe Work Load: <small>includes resident, support surface, head/footboards, and accessories</small>	600 lb	550 lb	450 lb	450 lb
Maximum Patient Weight: <small>EVENLY DISTRIBUTED</small>	550 lb	500 lb	400 lb	400 lb
Actuators:	4 DC	4 DC	4 DC	3 DC

## Options & Features:

Chair Position <small>(APS)</small> :	■	■	■	—
Advanced Positioning Pendant:	■	■	■	—
Pendant Attaches to Either Side of Bed:	■	■	■	—
Trendelenburg/Reverse Trendelenburg Positioning: <small>via APS staff control</small>	■	■	■	—
Slide-W-I-D-E™ Option Available: <small>between 35 - 42 inches</small>	■	■	—	—
Patented Super-Efficient Transmission: <small>less strain on the motors</small>	■	—	—	—
Hi-Low Actuators are Unloaded in Low Position:	■	■	■	—
Staff Control: <sup>1</sup>	■	■	■	—
Upgrade Hand Controls and/or Staff Controls Anytime: <small>interchangeable, mix and/or match standard and/or APS controls</small>	■	■	■	—
Splash Proof Electronics:	■	■	■	—
Intuitive Pedal Lock or 4 Locking Casters:	■	■	■	■ <sup>1</sup>
Bed Locks and Rolls at Any Height: <small>reducing care-giver back strain</small>	■	■	■	■
Silent Operating, Sealed Ball Bearings: <small>Maintenance-Free</small>	■	■	■	—
Taller Casters: <small>approximately 7 inches of clearance</small>	—	—	—	■
Head and Footboard Mounting:	Tall Tubular	Low Profile	Low Profile	Low Profile
Built-in Push Bar and Breakage Reduction:	■	—	—	—
Maximum Clearance: <small>over bed tables and lifts</small>	■	■	■	5.1 inch <sup>4</sup>
Can be Used as Transport Devices:	■	■	■	—
FDA Guideline Compliance: <small>for Bed Entrapment Risk Reductions</small>	■	■	■	—
Silver Solutions Powder Coating: <small>process provides antimicrobial surface protection</small>	■	■	■	—
Made in the USA by Skilled Craftsmen: <small>built to withstand trouble-free use for many years</small>	■	■	■	—

## Accessories:

Edema Foot Ratchet: <small>(elevates the foot end of the bed)</small>	■	■	■	—
Plastic Pan Deck:	■	■	■	—
Backup Battery: <small>(external or mounted onboard)</small>	■	■	■	—
Extension Kit: <small>(converts from 84 to 88 inch length)</small>	■	■	—	—
Pin-on Wide Kit:	■	■	■	—

## Warranty:

Main Frame & Welds	15 years	15 years	15 years	10 years
Control Box & Actuator Motors	5 years	5 years	4 years	2 years
Hand Control Pendant & Cabling	3 years	3 years	3 years	1 year
Staff Control & Cabling	3 years	3 years	3 years	—
Head & Footboards	1 year	1 year	1 year	1 year

Product Page Link

■ Standard	■ Available	■ Optional	— Not Applicable/Available

<sup>1</sup> anytime add-on or with purchase

<sup>3</sup> boards made with 3DL material - all styles available in 3DL - ask for details — custom finishes available upon request (additional cost may incur)

<sup>4</sup> fixed base

<sup>5</sup> "Hospital Bed System Dimensional and Assessment Guidance to Reduce Entrapment" details: [www.fda.gov](http://www.fda.gov)

<sup>6</sup> preserves the cured coating during the service life of the bed & inhibits the growth of organisms that can detract from the appearance and coating functionality. inhibits the growth of bacterial odors and resists stains caused by mildew and/or mold

GF Health Products, Inc. makes no independent representations or empirical statements about Silver Solutions. This information is reprinted from: TCI Powder Coatings, EPA Guidance for EPA Registration 40 CFR 152.25

GF Health Products, Inc. is an ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice. The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2017, GF Health Products, Inc. All Rights Reserved. Basic American Medical Products is a trademark of GF Health Products, Inc.

**grahamfield**

GF700011RevH23