



BARIATRIC PATIENT MOBILITY



REGENCY WHEELCHAIR

450 & 525 SERIES FIXED-BACK & RECLINING WHEELCHAIRS

Designed for optimal bariatric patient transport as well as hospital and clinical use, Regency 450 and 525 wheelchairs are built with a reinforced, heavy-duty frame and reinforced at all critical load-bearing points to support up to 450 lb (450 models) or 525 lb (525 models) users. With its multiple seat widths ranging from 24", 26", 28" to 30", with 20" seat depth, and seat-height settings of 15 1/2" (525 models) or 17 1/2" (450 and 525 models), this heavy-duty wheelchair enables institutions and caregivers to accommodate a large variety of patients.

The reclining wheelchair series features a pneumatic-assist reclining back to facilitate upright-to-recline positioning. Rear anti-tip brackets, elevating legrests, and removable headrest are standard.

Available with either a rear 20" x 1" wheel with chrome handrims and non-marking polyurethane tires (525 models) or 24" x 1" wheel with black handrims and non-marking polyurethane tires (450 and 525 models).

Features

- Robust, steel reinforced frame
- Non-marking polyurethane tires
- Comfortable, easy-to-clean black nylon upholstery
- Removable desk- and full-length arm styles with padded armrests
- Adjustable standard footrests or elevating legrest with durable black footplates
- 450 models: Weight Capacity 450 lb - EVENLY DISTRIBUTED
- 525 models: Weight Capacity 525 lb - EVENLY DISTRIBUTED
- Designed to meet the requirements of Medicare code: K0007 for fixed back models

Reclining Model Only

- Pneumatic-assist reclining back to facilitate full recline positioning
- Removable full support headrest standard
- Adjustable length elevating legrest
- Rear anti-tippers standard
- Designed to meet the requirements of Medicare code: K0007/E1226 (Recliner Back)
 - + E0966 (Head Rest Extension)
 - + E0971 (Anti-tipping device)
 - + K0195 (Elevating Legrests)

Warranty

- Five (5) year limited warranty



above: 525 reclining model with removable full arms (shown) and elevating legrests (shown). 20" x 1" rear wheel with chrome handrims option (shown) – desk arm and 24" x 1" rear wheel configurations available

center: 450 fixed-back model with removable desk arms (shown) and adjustable swingaway footrest (shown) – other configurations available

bottom left: 525 reclining model with removable desk arms (shown) and elevating legrests (shown). 24" x 1" rear wheel with black handrims option (shown) – full arm and 20" x 1" rear wheel with chrome handrims option configurations available



Regency Bariatric Fixed-Back Wheelchair

450 Series

| Item ¹ | Meets HCPCS | Overall Chair Width | | Overall Chair Depth | | Overall Chair Height | | Folded Chair Width | | Seat Width | | Seat Depth | | Seat Height factory setting | | Back Height | | Wheel ⁴ Diameter | | Caster ¹ Diameter | | Product ¹ Weight w/o front rigging | | Weight Capacity EVENLY DISTRIBUTED | | Detach. Full arms | Detach. Desk | Adjust. Swingaway Foot | Elevating Leg rests | Removable Head | Rear Anti-Tippers | Push Bar |
|-------------------|-------------|---------------------|-------|---------------------|--------|----------------------|-------|--------------------|-------|------------|-------|------------|-------|--------------------------------|-------|-------------|-------|-----------------------------|-------|------------------------------|-------|--|-------|---------------------------------------|--------|----------------------|--------------|------------------------|------------------------|----------------|-------------------|----------|
| | | inch | cm | inch | cm | inch | cm | inch | cm | inch | cm | inch | cm | inch | cm | inch | cm | inch | cm | inch | cm | lb | kg | lb | kg | | | | | | | |
| 450242020 | K0007 | 33 | 83.82 | 43.5 | 109.22 | 37 | 93.98 | 17 | 43.18 | 24 | 60.96 | 20 | 50.80 | 17.5 | 44.45 | 18.5 | 46.99 | 24 | 60.96 | 8 | 20.32 | 63 | 28.58 | 450 | 204.12 | — | ■ | ■ | — | — | ■ | ■ |
| 450262020 | K0007 | 35 | 88.90 | 43.5 | 109.22 | 37 | 93.98 | 17 | 43.18 | 26 | 66.04 | 20 | 50.80 | 17.5 | 44.45 | 18.5 | 46.99 | 24 | 60.96 | 8 | 20.32 | 63 | 28.58 | 450 | 204.12 | — | ■ | ■ | — | — | ■ | ■ |
| 450282020 | K0007 | 37 | 93.98 | 43.5 | 109.22 | 37 | 93.98 | 17 | 43.18 | 28 | 71.12 | 20 | 50.80 | 17.5 | 44.45 | 18.5 | 46.99 | 24 | 60.96 | 8 | 20.32 | 63 | 28.58 | 450 | 204.12 | — | ■ | ■ | — | — | ■ | ■ |
| 450302020 | K0007 | 39 | 99.06 | 43.5 | 109.22 | 37 | 93.98 | 17 | 43.18 | 30 | 76.20 | 20 | 50.80 | 17.5 | 44.45 | 18.5 | 46.99 | 24 | 60.96 | 8 | 20.32 | 63 | 28.58 | 450 | 204.12 | — | ■ | ■ | — | — | ■ | ■ |

■ Standard ■ Optional ■ Available ■ Not Available ■ Accessory — Does Not Apply 🔗 Active Hyperlink

¹ approximate product weight — as seat width increases, so does product weight, however slightly
³ adjustable height full arm accessory may be ordered separately at additional cost - item No. G10063

[450 Fixed Web Page](#)

Regency Bariatric Reclining Wheelchair

450 Series

| Item ¹ | Meets HCPCS | Overall Chair Width | | Overall Chair Depth | | Overall Chair Height | | Folded Chair Width | | Seat Width | | Seat Depth | | Seat Height factory setting | | Back Height | | Wheel ⁴ Diameter | | Caster ¹ Diameter | | Product ¹ Weight w/o front rigging | | Weight Capacity EVENLY DISTRIBUTED | | Detach. Full arms | Detach. Desk | Adjust. Swingaway Foot | Elevating Leg rests | Removable Head | Rear Anti-Tippers | Push Bar |
|-------------------|--------------------|---------------------|-------|---------------------|--------|----------------------|--------|--------------------|-------|------------|-------|------------|-------|--------------------------------|-------|-------------|-------|-----------------------------|-------|------------------------------|-------|--|-------|---------------------------------------|--------|----------------------|--------------|------------------------|------------------------|----------------|-------------------|----------|
| | | inch | cm | inch | cm | inch | cm | inch | cm | inch | cm | inch | cm | inch | cm | inch | cm | inch | cm | inch | cm | lb | kg | lb | kg | | | | | | | |
| 450242030R | multi ² | 33 | 83.82 | 52.5 | 133.35 | 52.5 | 133.35 | 17 | 43.18 | 24 | 60.96 | 20 | 50.80 | 17.5 | 44.45 | 36 | 91.44 | 24 | 60.96 | 8 | 20.32 | 94 | 42.64 | 450 | 204.12 | — | ■ | — | ■ | ■ | ■ | ■ |
| 450262030R | multi ² | 35 | 88.90 | 52.5 | 133.35 | 52.5 | 133.35 | 17 | 43.18 | 26 | 66.04 | 20 | 50.80 | 17.5 | 44.45 | 36 | 91.44 | 24 | 60.96 | 8 | 20.32 | 94 | 42.64 | 450 | 204.12 | — | ■ | — | ■ | ■ | ■ | ■ |
| 450282030R | multi ² | 37 | 93.98 | 52.5 | 133.35 | 52.5 | 133.35 | 17 | 43.18 | 28 | 71.12 | 20 | 50.80 | 17.5 | 44.45 | 36 | 91.44 | 24 | 60.96 | 8 | 20.32 | 94 | 42.64 | 450 | 204.12 | — | ■ | — | ■ | ■ | ■ | ■ |
| 450302030R | multi ² | 39 | 99.06 | 52.5 | 133.35 | 52.5 | 133.35 | 17 | 43.18 | 30 | 76.20 | 20 | 50.80 | 17.5 | 44.45 | 36 | 91.44 | 24 | 60.96 | 8 | 20.32 | 94 | 42.64 | 450 | 204.12 | — | ■ | — | ■ | ■ | ■ | ■ |

■ Standard ■ Optional ■ Available ■ Not Available ■ Accessory — Does Not Apply 🔗 Active Hyperlink

² designed to meet the requirements of Medicare code: K0007/E1226 (Recliner Back) + E0966 (Head Rest Extension) + E0971 (Anti-tipping device) + K095 (Elevating Legrests)
³ adjustable height full arm accessory may be ordered separately at additional cost - item No. G10063

[450 Reclining Web Page](#)

Regency Bariatric Fixed-Back Wheelchair
525 Series

| Item ¹ | Meets HCPCS | Overall Chair Width | | Overall Chair Depth | | Overall Chair Height | | Folded Chair Width | | Seat Width | | Seat Depth | | Seat Height factory setting | | Back Height | | Wheel ⁴ Diameter | | Caster ¹ Diameter | | Product ¹ Weight w/o front rigging | | Weight Capacity EVENLY DISTRIBUTED | | Detach. Full arms | Detach. Desk | Adjust. Swingaway Foot rests | Elevating Leg | Removable Head | Rear Anti-Tippers | Push Bar |
|-------------------|-------------|---------------------|-------|---------------------|--------|----------------------|--------|--------------------|-------|------------|-------|------------|-------|--------------------------------|-------|-------------|-------|-----------------------------|-------|------------------------------|-------|--|-------|---------------------------------------|--------|----------------------|--------------|---------------------------------|---------------|----------------|-------------------|----------|
| | | inch | cm | inch | cm | inch | cm | inch | cm | inch | cm | inch | cm | inch | cm | inch | cm | inch | cm | inch | cm | lb | kg | lb | kg | | | | | | | |
| 5252420520 | K0007 | 31 | 78.74 | 43.5 | 109.22 | 43.5 | 109.22 | 19 | 48.26 | 24 | 60.96 | 20 | 50.80 | 15.5 | 39.37 | 18.5 | 46.99 | 20 | 50.80 | 5 | 12.70 | 63 | 28.58 | 525 | 238.14 | — | ■ | ■ | — | — | ■ | ■ |
| 5252420540 | K0007 | 31 | 78.74 | 43.5 | 109.22 | 43.5 | 109.22 | 19 | 48.26 | 24 | 60.96 | 20 | 50.80 | 15.5 | 39.37 | 18.5 | 46.99 | 20 | 50.80 | 5 | 12.70 | 63 | 28.58 | 525 | 238.14 | ■ | — | ■ | — | — | ■ | ■ |
| 5252420720 | K0007 | 31 | 78.74 | 43.5 | 109.22 | 43.5 | 109.22 | 19 | 48.26 | 24 | 60.96 | 20 | 50.80 | 17.5 | 44.45 | 18.5 | 46.99 | 24 | 60.96 | 8 | 20.32 | 63 | 28.58 | 525 | 238.14 | — | ■ | ■ | — | — | ■ | ■ |
| 5252420740 | K0007 | 31 | 78.74 | 43.5 | 109.22 | 43.5 | 109.22 | 19 | 48.26 | 24 | 60.96 | 20 | 50.80 | 17.5 | 44.45 | 18.5 | 46.99 | 24 | 60.96 | 8 | 20.32 | 63 | 28.58 | 525 | 238.14 | ■ | — | ■ | — | — | ■ | ■ |
| 5252620520 | K0007 | 33 | 83.82 | 43.5 | 109.22 | 43.5 | 109.22 | 19 | 48.26 | 26 | 66.04 | 20 | 50.80 | 15.5 | 39.37 | 18.5 | 46.99 | 20 | 50.80 | 5 | 12.70 | 63 | 28.58 | 525 | 238.14 | — | ■ | ■ | — | — | ■ | ■ |
| 5252620540 | K0007 | 33 | 83.82 | 43.5 | 109.22 | 43.5 | 109.22 | 19 | 48.26 | 26 | 66.04 | 20 | 50.80 | 15.5 | 39.37 | 18.5 | 46.99 | 20 | 50.80 | 5 | 12.70 | 63 | 28.58 | 525 | 238.14 | ■ | — | ■ | — | — | ■ | ■ |
| 5252620720 | K0007 | 33 | 83.82 | 43.5 | 109.22 | 43.5 | 109.22 | 19 | 48.26 | 26 | 66.04 | 20 | 50.80 | 17.5 | 44.45 | 18.5 | 46.99 | 24 | 60.96 | 8 | 20.32 | 63 | 28.58 | 525 | 238.14 | — | ■ | ■ | — | — | ■ | ■ |
| 5252620740 | K0007 | 33 | 83.82 | 43.5 | 109.22 | 43.5 | 109.22 | 19 | 48.26 | 26 | 66.04 | 20 | 50.80 | 17.5 | 44.45 | 18.5 | 46.99 | 24 | 60.96 | 8 | 20.32 | 63 | 28.58 | 525 | 238.14 | ■ | — | ■ | — | — | ■ | ■ |
| 5252820520 | K0007 | 35 | 88.90 | 43.5 | 109.22 | 43.5 | 109.22 | 19 | 48.26 | 28 | 71.12 | 20 | 50.80 | 15.5 | 39.37 | 18.5 | 46.99 | 20 | 50.80 | 5 | 12.70 | 63 | 28.58 | 525 | 238.14 | — | ■ | ■ | — | — | ■ | ■ |
| 5252820540 | K0007 | 35 | 88.90 | 43.5 | 109.22 | 43.5 | 109.22 | 19 | 48.26 | 28 | 71.12 | 20 | 50.80 | 15.5 | 39.37 | 18.5 | 46.99 | 20 | 50.80 | 5 | 12.70 | 63 | 28.58 | 525 | 238.14 | ■ | — | ■ | — | — | ■ | ■ |
| 5252820720 | K0007 | 35 | 88.90 | 43.5 | 109.22 | 43.5 | 109.22 | 19 | 48.26 | 28 | 71.12 | 20 | 50.80 | 17.5 | 44.45 | 18.5 | 46.99 | 24 | 60.96 | 8 | 20.32 | 63 | 28.58 | 525 | 238.14 | — | ■ | ■ | — | — | ■ | ■ |
| 5252820740 | K0007 | 35 | 88.90 | 43.5 | 109.22 | 43.5 | 109.22 | 19 | 48.26 | 28 | 71.12 | 20 | 50.80 | 17.5 | 44.45 | 18.5 | 46.99 | 24 | 60.96 | 8 | 20.32 | 63 | 28.58 | 525 | 238.14 | ■ | — | ■ | — | — | ■ | ■ |
| 5253020520 | K0007 | 37 | 93.98 | 43.5 | 109.22 | 43.5 | 109.22 | 19 | 48.26 | 30 | 76.20 | 20 | 50.80 | 15.5 | 39.37 | 18.5 | 46.99 | 20 | 50.80 | 5 | 12.70 | 63 | 28.58 | 525 | 238.14 | — | ■ | ■ | — | — | ■ | ■ |
| 5253020540 | K0007 | 37 | 93.98 | 43.5 | 109.22 | 43.5 | 109.22 | 19 | 48.26 | 30 | 76.20 | 20 | 50.80 | 15.5 | 39.37 | 18.5 | 46.99 | 20 | 50.80 | 5 | 12.70 | 63 | 28.58 | 525 | 238.14 | ■ | — | ■ | — | — | ■ | ■ |
| 5253020720 | K0007 | 37 | 93.98 | 43.5 | 109.22 | 43.5 | 109.22 | 19 | 48.26 | 30 | 76.20 | 20 | 50.80 | 17.5 | 44.45 | 18.5 | 46.99 | 24 | 60.96 | 8 | 20.32 | 63 | 28.58 | 525 | 238.14 | — | ■ | ■ | — | — | ■ | ■ |
| 5253020740 | K0007 | 37 | 93.98 | 43.5 | 109.22 | 43.5 | 109.22 | 19 | 48.26 | 30 | 76.20 | 20 | 50.80 | 17.5 | 44.45 | 18.5 | 46.99 | 24 | 60.96 | 8 | 20.32 | 63 | 28.58 | 525 | 238.14 | ■ | — | ■ | — | — | ■ | ■ |

■ Standard ■ Optional ■ Available ■ Not Available ■ Accessory — Does Not Apply 🔗 Active Hyperlink

¹ approximate product weight — as seat width increases, so does product weight, however slightly
³ adjustable height full arm accessory may be ordered separately at additional cost - Item No. G10063
⁴ chairs ending in 520, 530, 540 and 550 all have 20" rear wheels and 5" plastic front casters.

[525 Web Page](#)

Regency Bariatric Reclining Wheelchair
525 Series

| Item* | Meets HCPCS | Overall Chair Width | | Overall Chair Depth | | Overall Chair Height | | Folded Chair Width | | Seat Width | | Seat Depth | | Seat Height factory setting | | Back Height | | Wheel ⁴ Diameter | | Caster ¹ Diameter | | Product ¹ Weight w/o front rigging | | Weight Capacity EVENLY DISTRIBUTED | | Detach. Full arms | Detach. Desk | Adjust. Swingaway Foot rests | Elevating Leg | Removable Head | Rear Anti-Tippers | Push Bar | |
|-------------|--------------------|---------------------|--------|---------------------|--------|----------------------|--------|--------------------|-------|------------|-------|------------|-------|-----------------------------|-------|-------------|------|-----------------------------|-------|------------------------------|-------|---|-------|------------------------------------|--------|-------------------|--------------|------------------------------|---------------|----------------|-------------------|----------|---|
| | | inch | cm | inch | cm | inch | cm | inch | cm | inch | cm | inch | cm | inch | cm | inch | cm | inch | cm | inch | cm | lb | kg | lb | kg | | | | | | | | |
| 5252420530R | multi ² | 34 | 86.36 | 52.5 | 133.35 | 52.5 | 133.35 | 19 | 48.26 | 24 | 60.96 | 20 | 50.80 | 15.5 | 39.37 | 36 | 91.4 | 20 | 50.80 | 8 | 20.32 | 94 | 42.64 | 525 | 238.14 | — | ■ | — | ■ | ■ | ■ | ■ | ■ |
| 5252420550R | multi ² | 34 | 86.36 | 52.5 | 133.35 | 52.5 | 133.35 | 19 | 48.26 | 24 | 60.96 | 20 | 50.80 | 15.5 | 39.37 | 36 | 91.4 | 20 | 50.80 | 5 | 12.70 | 94 | 42.64 | 525 | 238.14 | ■ | — | — | ■ | ■ | ■ | ■ | ■ |
| 5252420730R | multi ² | 34 | 86.36 | 52.5 | 133.35 | 52.5 | 133.35 | 19 | 48.26 | 24 | 60.96 | 20 | 50.80 | 17.5 | 44.45 | 36 | 91.4 | 24 | 60.96 | 8 | 20.32 | 94 | 42.64 | 525 | 238.14 | — | ■ | — | ■ | ■ | ■ | ■ | ■ |
| 5252420750R | multi ² | 34 | 86.36 | 52.5 | 133.35 | 52.5 | 133.35 | 19 | 48.26 | 24 | 60.96 | 20 | 50.80 | 17.5 | 44.45 | 36 | 91.4 | 24 | 60.96 | 8 | 20.32 | 94 | 42.64 | 525 | 238.14 | ■ | — | — | ■ | ■ | ■ | ■ | ■ |
| 5252620530R | multi ² | 36 | 91.44 | 52.5 | 133.35 | 52.5 | 133.35 | 19 | 48.26 | 26 | 66.04 | 20 | 50.80 | 15.5 | 39.37 | 36 | 91.4 | 20 | 50.80 | 5 | 12.70 | 94 | 42.64 | 525 | 238.14 | — | ■ | — | ■ | ■ | ■ | ■ | ■ |
| 5252620550R | multi ² | 36 | 91.44 | 52.5 | 133.35 | 52.5 | 133.35 | 19 | 48.26 | 26 | 66.04 | 20 | 50.80 | 15.5 | 39.37 | 36 | 91.4 | 20 | 50.80 | 5 | 12.70 | 94 | 42.64 | 525 | 238.14 | ■ | — | — | ■ | ■ | ■ | ■ | ■ |
| 5252620730R | multi ² | 36 | 91.44 | 52.5 | 133.35 | 52.5 | 133.35 | 19 | 48.26 | 26 | 66.04 | 20 | 50.80 | 17.5 | 44.45 | 36 | 91.4 | 24 | 60.96 | 8 | 20.32 | 94 | 42.64 | 525 | 238.14 | — | ■ | — | ■ | ■ | ■ | ■ | ■ |
| 5252620750R | multi ² | 36 | 91.44 | 52.5 | 133.35 | 52.5 | 133.35 | 19 | 48.26 | 26 | 66.04 | 20 | 50.80 | 17.5 | 44.45 | 36 | 91.4 | 24 | 60.96 | 8 | 20.32 | 94 | 42.64 | 525 | 238.14 | ■ | — | — | ■ | ■ | ■ | ■ | ■ |
| 5252820530R | multi ² | 38 | 96.52 | 52.5 | 133.35 | 52.5 | 133.35 | 19 | 48.26 | 28 | 71.12 | 20 | 50.80 | 15.5 | 39.37 | 36 | 91.4 | 20 | 50.80 | 5 | 12.70 | 94 | 42.64 | 525 | 238.14 | — | ■ | — | ■ | ■ | ■ | ■ | ■ |
| 5252820550R | multi ² | 38 | 96.52 | 52.5 | 133.35 | 52.5 | 133.35 | 19 | 48.26 | 28 | 71.12 | 20 | 50.80 | 15.5 | 39.37 | 36 | 91.4 | 20 | 50.80 | 5 | 12.70 | 94 | 42.64 | 525 | 238.14 | ■ | — | — | ■ | ■ | ■ | ■ | ■ |
| 5252820730R | multi ² | 38 | 96.52 | 52.5 | 133.35 | 52.5 | 133.35 | 19 | 48.26 | 28 | 71.12 | 20 | 50.80 | 17.5 | 44.45 | 36 | 91.4 | 24 | 60.96 | 8 | 20.32 | 94 | 42.64 | 525 | 238.14 | — | ■ | — | ■ | ■ | ■ | ■ | ■ |
| 5252820750R | multi ² | 38 | 96.52 | 52.5 | 133.35 | 52.5 | 133.35 | 19 | 48.26 | 28 | 71.12 | 20 | 50.80 | 17.5 | 44.45 | 36 | 91.4 | 24 | 60.96 | 8 | 20.32 | 94 | 42.64 | 525 | 238.14 | ■ | — | — | ■ | ■ | ■ | ■ | ■ |
| 5253020530R | multi ² | 40 | 101.60 | 52.5 | 133.35 | 52.5 | 133.35 | 19 | 48.26 | 30 | 76.20 | 22 | 55.88 | 15.5 | 39.37 | 36 | 91.4 | 20 | 50.80 | 5 | 12.70 | 94 | 42.64 | 525 | 238.14 | — | ■ | — | ■ | ■ | ■ | ■ | ■ |
| 5253020550R | multi ² | 40 | 101.60 | 52.5 | 133.35 | 52.5 | 133.35 | 19 | 48.26 | 30 | 76.20 | 22 | 55.88 | 15.5 | 39.37 | 36 | 91.4 | 20 | 50.80 | 5 | 12.70 | 94 | 42.64 | 525 | 238.14 | ■ | — | — | ■ | ■ | ■ | ■ | ■ |
| 5253020730R | multi ² | 40 | 101.60 | 52.5 | 133.35 | 52.5 | 133.35 | 19 | 48.26 | 30 | 76.20 | 22 | 55.88 | 17.5 | 44.45 | 36 | 91.4 | 24 | 60.96 | 8 | 20.32 | 94 | 42.64 | 525 | 238.14 | — | ■ | — | ■ | ■ | ■ | ■ | ■ |
| 5253020750R | multi ² | 40 | 101.60 | 52.5 | 133.35 | 52.5 | 133.35 | 19 | 48.26 | 30 | 76.20 | 22 | 55.88 | 17.5 | 44.45 | 36 | 91.4 | 24 | 60.96 | 8 | 20.32 | 94 | 42.64 | 525 | 238.14 | ■ | — | — | ■ | ■ | ■ | ■ | ■ |

■ Standard ■ Optional ■ Available ■ Not Available ■ Accessory — Does Not Apply 🔗 Active Hyperlink

¹ approximate product weight — as seat width increases, so does product weight, however slightly
² designed to meet the requirements of Medicare code: K0007/E1226 (Recliner Back) + E0966 (Head Rest Extension) + E0971 (Anti-tipping device) + K095 (Elevating Legrests)
³ adjustable height full arm accessory may be ordered separately at additional cost - item No. G10063
⁴ chairs ending in 520, 530, 540 and 550 all have 20" rear wheels and 5" plastic front casters.

🔗 525 Web Page

GF Health Products, Inc. is an ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice. The most current and complete product information can be found on our website.

www.gendroninc.com | www.grahamfield.com



© 2021, GF Health Products, Inc. All Rights Reserved. Gendron is a registered trademark of GF Health Products, Inc.

graham field

GFz100014 Rev C24