

## Features:

	LS200	LS300	LS400
True Low Air Loss Function:	■ 4.5L / min.	■ 6L / min.	■ 6L / min.
2:1 Alternating Function:	■	■	■
Half Body: <small>always supported at any one time</small>	■	■	■
Cell on Cell Therapy:	■	■	■
Static Function :	■	■	■
Low Pressure Alarm:	audio/visual	audio/visual	audio/visual
Inflation Time:	approx 45 min	Fast Inflate: approx 20 min <sup>3</sup>	Fast Inflate: approx 20 min <sup>3</sup>
Auto Firm: <sup>1</sup> <small>for Nursing Procedures</small>	■	■	■
Fowler Function: <sup>2</sup> <small>helps from bottoming out in Fowler position</small>	■	■	■
Soft/Firm Levels:	variable knob	5 digital	8 digital
Alternating Function Times: <sup>4</sup>	12 min	10, 15, 20 min	5, 10, 15, 20 min

## Mattress:

	inch		cm +		inch		cm +		inch		cm +		
	Width:	36	91.44	36	91.44	36	91.44	36	91.44	36	91.44	36	91.44
Length:	80	203.20	80	203.20	80	203.20	80	203.20	80	203.20	80	203.20	
Thickness:	8	20.32	8	20.32	8	20.32	8	20.32	8	20.32	8	20.32	
	lb		kg +		lb		kg +		lb		kg +		
	Weight:	20	9.07	16.3	7.39	17.7	8.03	17.7	8.03	17.7	8.03	17.7	8.03
<b>Cover:</b>	Quilted:	■	■	■	■	■	■	■	■	■	■	■	
	Zipper:	■	■	■	■	■	■	■	■	■	■	■	
	Cell Qty:	18	18	18	18	18	18	18	18	18	18	18	
<b>Material:</b>	Top:	nylon w/pu laminated		nylon w/pu laminated		nylon w/pu laminated		nylon w/pu laminated		nylon w/pu laminated		nylon w/pu laminated	
	Bottom:	polyester laminated with PVC		polyester laminated with PVC		polyester laminated with PVC		polyester laminated with PVC		polyester laminated with PVC		polyester laminated with PVC	

## Pump:

	inch		cm +		inch		cm +		inch		cm +	
	Width:	10.2	25.91	10.2	25.91	7.7	19.56	7.7	19.56	7.7	19.56	7.7
Length:	4.6	11.68	4.1	10.41	3	7.62	3	7.62	3	7.62	3	7.62
Height:	6.4	16.26	6.4	16.26	4.8	12.19	4.8	12.19	4.8	12.19	4.8	12.19
	lb		kg +		lb		kg +		lb		kg +	
	Weight:	4.9	2.223	4.9	2.223	7.7	3.493	7.7	3.493	7.7	3.493	7.7
Cycle Time:	12 min		10, 15, 20 min		5, 10, 15, 20 min		5, 10, 15, 20 min		5, 10, 15, 20 min		5, 10, 15, 20 min	
Flow Rate: <small>direct from pump</small>	■ 4.5L / min. compressor		■ 6L / min. compressor		■ 6L / min. compressor		■ 6L / min. compressor		■ 6L / min. compressor		■ 6L / min. compressor	
Pressure:	16~32 (± 6) mmHg		16~32 (± 6) mmHg		16~50 (± 6) mmHg		16~50 (± 6) mmHg		16~50 (± 6) mmHg		16~50 (± 6) mmHg	
Rated Voltage & Frequency:	AC 110-120V 60 Hz		AC 110-120V 60 Hz		AC 110-120V 60 Hz		AC 110-120V 60 Hz		AC 110-120V 60 Hz		AC 110-120V 60 Hz	
Fuse Rating:	T1AL / 250V		T1AL / 250V		T1AL / 250V		T1AL / 250V		T2AL / 250V		T2AL / 250V	
Maximum Current:	0.1A		0.1A		0.1A		0.1A		0.5A		0.5A	
AC Power Cable:	12.5 ft – non-shielding		12.5 ft – non-shielding		12.5 ft – non-shielding		12.5 ft – non-shielding		15 ft – non-shielding		15 ft – non-shielding	
Expected Service Life:	3 years		3 years		3 years		3 years		3 years		3 years	

<sup>1</sup> allows patient to be turned faster if needed; easier sling placement  
<sup>2</sup> increases mattress firmness by 2 levels to help prevent bottoming out in chair position  
<sup>3</sup> helpful for setup and in cases where nurses needs to perform certain procedures  
<sup>4</sup> more options provide healthcare professionals to better control bed ulcers and air flow  
 † standard rounding method used to show metric values as whole numbers

■ Standard ■ Not Available