


## LIFT SLINGS

# 2-POINT SLINGS

### Features

- One piece canvas sling for use with 2-point spreader
- Provides support from shoulders to knees
- Used for transferring, weighing or toileting
- For use with chain set, item GF133-S-C (not included with sling, sold separately)
- Designed to meet the requirements of HCPCS Code: E0621
- Six (6) month limited warranty



		body Width		body Length		commode opening Width		commode opening Length		Weight Capacity		UOM
		inch	cm	inch	cm	inch	cm	inch	cm	lb	Kg	
GF112-C-LC	Without Commode Opening	26	66.04	41.5	105.41					220	99.8	1 ea
GF113-C-LC	With Commode Opening	26	66.04	41.5	105.41	6.75	17.15	10.625	26.99	220	99.8	1 ea



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2007, GF Health Products, Inc. All Rights Reserved.  
Lumex is a registered trademark of GF Health Products, Inc.