




LIFT SLINGS

DELUXE
SIT-TO-STAND PADDED SLINGS

Features

- Cushioned support panels under the arm and lumbar areas enhance patient comfort
- Polyester with nylon stitching and webbing
- The unique cross-over strap and non-slip inner pad keep the sling in place so that it does not ride up
- Optional stand assist buttocks strap is available for initial lift assistance
- For use with LF2020 STS lift , DSLSA4 for use with LF2090 STS lift
- Six (6) month limited warranty

	Size	Torso Girth		Weight Capacity EVENLY DISTRIBUTED		UOM
		inch	cm	lb	Kg	
DSLSA1	Small	26 – 43	11.8 – 19.5	400	181.4	1 ea
DSLSA2	Medium	29 – 48	13.2 – 21.8	400	181.4	1 ea
DSLSA3	Large	36 – 59	16.3 – 26.8	400	181.4	1 ea
DSLSA4	XX-Large	50 – 79	22.7 – 35.8	600	272.2	1 ea



← Torso Girth →



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2007, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.