




LIFT SLINGS

PADDED TOILETING SLING

Features

- Provides an open area from the middle of the back to the middle of the thigh
- Polyester with nylon stitching and webbing
- Ideal for changing incontinence pads or toileting
- Designed to meet the requirements of HCPCS Code:E0621
- Maximum Weight Capacity 400 lb - EVENLY DISTRIBUTED
- Six (6) month limited warranty



Size	body Width		body Length		total Width		total Length		Weight Capacity EVENLY DISTRIBUTED		UOM	
	inch	cm	inch	cm	inch	cm	inch	cm	lb	Kg		
TP332	Medium	33.25	84.46	11.5	29.21	50.75	128.91	56.25	142.88	400	181.4	1 ea
TP333	Large	36.5	92.71	14	35.56	56.5	143.51	65.5	166.37	400	181.4	1 ea



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2010, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.

graham
field

GF1000076RevB13